
Varillas magnéticas de agitación Spinbar (continuación)

Varillas de forma ovalada






- Adaptado a los recipientes y matraces de fondo redondo

An x Ø (mm)	Para matraces	Ref.	€
10 x 5	-	045982	NC -
15 x 6	-	045983	NC -
19 x 9,5	De 10 ml a 100 ml	045984	NC -
25 x 12	De 100 ml a 200 ml	045985	NC -
30 x 16	De 300 ml a 500 ml	045986	NC -
35 x 16	De 300 ml a 500 ml	045987	NC -
40 x 13	De 500 ml a 1000 ml	046097	NC -
40 x 20	De 500 ml a 1000 ml	045988	NC -
50 x 20	De 2000 ml a 3000 ml	045989	NC -
64 x 20	De 3000 ml a 5000 ml	045990	NC -
70 x 20	De 12 l a 22 l	045991	NC -

Varillas triangulares





- Especialmente adapatadas para disolver sólidos y sedimentos
- La fricción de la base provoca una acción de rascado en el fondo


An x Ø (mm)	Ref.	€	An x Ø (mm)	Ref.	€
12 x 6	095882	NC -	50 x 12	085887	NC -
20 x 8	095883	NC -	55 x 14	085888	NC -
25 x 8	095884	NC -	80 x 18	085889	NC -
25 x 14	046087	NC -	110 x 36	085890	NC -
35 x 9	095885	NC -	136 x 36	085891	NC -
40 x 14	085886	NC -	-	-	-

Varillas cuadradas

- Fuertes turbulencias, sobre todo a baja velocidad

An x Ø (mm)	Ref.	€
12 x 4	095878	NC -
25 x 5,5	095879	NC -
35 x 6	095880	NC -
50 x 7,5	095881	NC -







Varillas disco

- Para utilizar en tubos de ensayo


An x Ø (mm)	Ref.	€
6 x 9	098889	NC -
6 x 10	095890	NC -
10 x 20	098891	NC -
12 x 30	095892	NC -





Varillas perfiladas






- Acción similar a las ovaladas o pivotantes

An x Ø (mm)	Ref.	€	An x Ø (mm)	Ref.	€
10 x 4	095931	NC -	40 x 8	095937	NC -
15 x 5	095932	NC -	45 x 8	095938	NC -
20 x 7	095933	NC -	50 x 8	095939	NC -
25 x 8	095934	NC -	60 x 8	095941	NC -
30 x 8	095935	NC -	70 x 8	095942	NC -
35 x 8	095936	NC -	80 x 10	095943	NC -

Varillas con pivote extraíble





- Puede utilizarse con o sin pivote
- Perfectamente cilíndricas y lisas


An x Ø (mm)	Ref.	€	An x Ø (mm)	Ref.	€
12 x 8	095919	NC -	51 x 8	095925	NC -
25 x 8	095920	NC -	51 x 10	095926	NC -
32 x 8	095921	NC -	75 x 12	095927	NC -
38 x 8	095922	NC -	102 x 16	095928	NC -
38 x 10	095923	NC -	127 x 16	095929	NC -
42 x 10	095924	NC -	150 x 19	095930	NC -
45 x 8	095975	NC -	-	-	-

Varillas cilíndricas de vidrio

- Evita el contacto con el PTFE

An x Ø (mm)	Ref.	€
6 x 5	095910	NC -
12 x 5	095911	NC -
25 x 6	095912	NC -
45 x 8	095913	NC -
60 x 8	095914	NC -






Varillas con forma de cruz

- Muy estable, provoca fuertes turbulencias

An x Ø (mm)	Ref.	€
10 x 5	095951	NC -
20 x 8	095952	NC -
25 x 9	095953	NC -
30 x 10	095954	NC -
38 x 11	095955	NC -
50 x 15	095956	NC -
60 x 20	095957	NC -





Precios indicativos. No incluyen IVA. Visiten nuestras WEB