
CENTRIFUGE TUBES / CULTURE TUBES WITH DUAL POSITION CAP

Nunc™ large volume centrifuge tube



Volume (ml)	250
Ø (mm)	59.7
Height (mm)	144.2
Material	Polypropylene; polyethylene cap
Sterile	Yes
Bottom	Conical
Autoclavable	Yes
Max centrifugal speed (g)	10000
Graduation	No
Marking area	No
Specifications	Free of DNase, RNase and endotoxin; USP Class VI compliant
Units/carton	10 bags of 40 tubes
Blue cap Cat. No.	**056886**
€ Excl. VAT/carton	**NC -**

Falcon®Culture tubes with dual position cap



Volume (ml)	5							14		
Ø ext. (mm)	12.4	12.4	12.4	12.4	12.4	12.4	12.4	17.6	17.6	17.6
Ø int. (mm)	10.72	10.72	10.72	10.72	10.62	10.62	10.62	15.77	15.77	15.77
Height without cap (mm)	75.54	75.54	75.54	75.54	73.79	73.79	73.79	95.58	94.82	94.82
Height with cap (mm)	76.22	76.22	-	-	74.54	-	-	96.09	96.09	-
Material	Polystyrene; LDPE cap	Polystyrene; LDPE cap	Polystyrene	Polystyrene	Polypropylene; LDPE cap	Polypropylene	Polypropylene	Polystyrene; LDPE cap	Polypropylene; LDPE cap	Polypropylene
Sterile	Yes	Yes	Yes	No	Yes	Yes	No	Yes	Yes	Yes
Bottom	Rounded	Rounded	Rounded	Rounded	Rounded	Rounded	Rounded	Rounded	Rounded	Rounded
Caps	Yes	Yes	No	No	Yes	No	No	Yes	Yes	No
Type of cap	Snap; Dual position	Snap; Dual position with 35 µm sieve	-	-	Snap; Dual position	-	-	Snap; Dual position	Snap; Dual position	-
Autoclavable	No	No	No	No	No	No	No	No	No	No
Max centrifugal speed (g)	1400	1400	1400	1400	3000	3000	3000	1400	3000	3000
Graduation	No	No	No	No	No	No	No	Yes; moulded	Yes; screen printed	Yes; screen printed
Marking area	No	No	No	No	No	No	No	No	Yes	Yes
Colour	Transparent	Transparent	Transparent	Transparent	Translucent	Translucent	Translucent	Transparent	Translucent	Translucent
Operating temperature	+4 °C to +60 °C	+4 °C to +60 °C	+4 °C to +60 °C	+4 °C to +60 °C	-80 °C to +121 °C	-80 °C to +121 °C	-80 °C to +121 °C	+4 °C to +60 °C	-80 °C to +121 °C	-80 °C to +121 °C
Specifications	Free of DNase, RNase and endotoxin <0,1EU/ml; Hydrophobic and biologically inert surface			Free of DNase, RNase ; Hydrophobic and biologically inert surface	Free of DNase, RNase and endotoxin <0,1EU/ml; Hydrophobic and biologically inert surface		Free of DNase, RNase ; Hydrophobic and biologically inert surface	Free of DNase, RNase and endotoxin <0,1EU/ml; Hydrophobic and biologically inert surface		
Cat. No.	**352003**	-	-	-	-	-	-	**352001**	**352006**	-
Units/carton	500 individually wrapped tubes									
€ Excl. VAT/carton	**NC -**	-	-	-	-	-	-	**NC -**	**NC -**	-
Cat. No.	**352054**	-	**352052**	-	-	**352053**	-	**352051**	-	**352018**
Units/carton	8 bags of 125 tubes									
€ Excl. VAT/carton	**NC -**	-	**NC -**	-	-	**NC -**	-	**NC -**	-	**NC -**
Cat. No.	**352058**	**352235**	-	-	**352063**	-	-	**352057**	**352059**	-
Units/carton	20 blister of 25 tubes									
€ Excl. VAT/carton	**NC -**	**NC -**	-	-	**NC -**	-	-	**NC -**	**NC -**	-
Cat. No.	-	-	-	**352008**	-	-	**352002**	-	-	-
Units/carton	Bag of 1000 tubes									
€ Excl. VAT/carton	-	-	-	**NC -**	-	-	**NC -**	-	-	-

Dual position snap plug

Cat. No.	**352032**
Units/carton	2000
€ Excl. VAT/carton	**NC -**

