Midigel XL horizontal tank



Cat. No.	Description	€ Excl. VAT
049416	Midigel XL Tank complete	NC -
049427	UV transparent Tray 15 x 17 cm (1 to 6 combs)	NC -
049428	UV transparent Tray 15 x 12 cm (1 to 4 combs)	NC -
049429	UV transparent Tray 15 x 7 cm (1 to 2 combs)	NC -
049003	Set of cords	NC -
049004	Separate gel casting system with lwith adjustable level and feet	NC -

Midigel XL tank packs for 96 samples

Pack content	Midigel XL tank + gel holder 150 x 170 mm + gel holder 150 x 120 mm + removable walls for direct gel pouring + 6 combs 16 teeth, thickness 1 mm for multichannel pipettes + safety lid + set of cords, user manual	
Generator included	No	PS 304XL (300 V - 400 mA - 100 W)
Cat.No.	**049903B**	**049902B**
€ Excl. VAT	**NC -**	**NC -**

Gel size : 15 x 17 cm.

- **Buffer volume : 1200 ml**
- **Maximum number of samples : 192 (with 6 combs)**
- **Dimensions : 32 x 23 x 9 cm**
- **For standard size gels of 17, 12 and 7 cm in length depending on the rack used**
- **For PCR and RFLP screening**
- **For analytical or preparative applications on nucleic acids, agarose gel restriction fragments**
- **Gel tray 15 x 17 cm with 6 notches for comb insertion**
- **Gel tray 15 x 12 cm with 4 notches**
- **Gel tray 15 x 7 cm with 2 notches**
- **The holders are UV-transparent and equipped with sliding walls for gel casting**
- **16 and 32 tooth combs are compatible with multi-channel pipettes**
- **Preparatory combs available**
- **Adjustable height feet and built-in spirit level**
- **Red strips deposits visualisation**
- **The tank is delivered complete**

Suggested generators

- PS 304 : 300 Volts - 400 mA - 120 Watts
- PS 305 : 300 Volts - 500 mA - 150 Watts with Timer and Gel Saver
- PS 503 : 500 Volts - 300 mA - 150 Watts with Timer and Gel Saver

The tank is supplied with :

- Removable UV transparent gel tray 15 x 17 cm;
- 2 removable walls for direct casting of gels;
- 1 comb, 2 mm thickness , 27 teeth, volume 32 µl ;
- safety cover ;
- set of cords ;
- user's manual.

Maxigel ECO 2 horizontal tank



Cat. No.	Description	€ Excl. VAT
049039	Maxigel ECO 2	NC -
049040	UV-transparent tray 20 x 24 cm	NC -
049041	UV transparent Tray 20 x 12 cm	NC -
049003	Set of cords	NC -
049004	Separate gel casting system with lwith adjustable level and feet	NC -



Find all compatible combs on www.dulis.nl

Gel size : 20 x 24 cm or 20 x 12 cm.

- **Buffer volume : 1800 ml**
- **Max. number of samples : 352 (with 8 combs 049053 or 049461 compatible with multichannel pipettes)**
- **Dimensions : 28 x 37 x 9 cm**
- **For DNA or RNA analysis**
- **For rapid screening after PCR of 96 samples on a single gel**
- **For restriction fragment screening in agarose gel**
- **Gel tray 20 x 24 cm with eight notches for inserting the combs**
- **Gel tray 20 x 12 cm with 3 notches**
- **44-tooth combs compatible with multi-channel pipettes**
- **Preparatory combs available**
- **Adjustable height feet and built-in spirit level**
- **Red strips deposits visualisation**
- **The tank is delivered complete**

Suggested generators

- PS 305 : 300 Volts - 500 mA - 150 Watts with Timer and Gel Saver
- PS 503 : 500 Volts - 300 mA - 150 Watts with Timer and Gel Saver
- PS 608 : 600 Volts - 800 mA - 300 Watts with Timer and 3 displays

The tank is supplied with :

- Removable UV transparent gel tray 20 x 24 cm;
- 2 removable walls for direct casting of gels;
- 1 comb, 1 mm thickness, 44/21 teeth - volume 14/24 µl (Cat.No. 049053)
- safety cover ;
- set of cords ;
- accessories for buffer recirculation ;
- user's manual.