Pipet-Aid®









Pipette controller and reagent dispenser for glass and plastic pipettes from 1 to 100 ml. The Pipet-Aid is the most widely used pipette controler in the world: its reliability, robustness and ease of maintenance make it the most renowned instrument in its category.

4 models form a complete family:

- Pipet-Aid XP2: ergonomic design; comfortable to hold; charging while handling, thanks to the new charger; ultra-quiet pump; supplied with 4 additional filters, holder and charger;
- Pipet-Aid: is the basic instrument, light, autonomous, wireless;
- Pipet-Aid XP: faster, more compact, more functions, delivered with support ;
- Pipet-Aid XP Gravity: identical to Pipet-Aid XP with an additional gravity function (gravity dispensing, without pushing). Gravity system, useful in cell culture to avoid peeling off a cell mat.

	Pipet-Aid XP2	Pipet-Aid	Pipet-Aid XP	Pipet-Aid XP Gravity
Rechargeable battery 9 V	Yes	Yes	Yes	Yes
Filtering nosepiece for culture	Yes	Yes	Yes	Yes
Speed	3 speeds*	1 speed (normal)	3 speeds*	3 speeds* + gravity
Pipettes accepted	1 to 100 ml	1 to 100 ml	1 to 100 ml	1 to 100 ml
Cat. No	**075212**	**075001**	**075002**	**075003**
€ Excl.VAT of Pipet-Aid (complete with charger)	**NC -**	**NC -**	**NC -**	**NC -**

* (normal, slow, fast)

Spare nose

Colours	Black	Rede	Green	Blue	Orange	Violet	Yellow	Assorted	Nosepiece for XP2
Cat. No.	**075004**	**075005**	**075006**	**075007**	**075008**	**075009**	**075109**	**075010**	**075213**
Units/box	1	2	2	2	2	2	2	6	1
€ Excl.VAT/box	**NC -**	**NC -**	**NC -**	**NC -**	**NC -**	**NC -**	**NC -**	**NC -**	**NC -**

Spare parts

Cat. No.	Description	€ Excl. VAT
	Parts common to all 4 Pipet-Aid models	
075011	Hydrophobic filter 0.8 µm (5 units)	**NC -**
075233	Hydrophobic sterile filter 0,8 µm (5 units)	**NC -**
075015	Silicone end pipette holder (2 units)	**NC -**
	Common parts to Wireless Pipet-Aid and Pipet-Aid Hood Mate	
075012	Filter holding gasket (2 units)	**NC -**
075013	Battery 9 V	**NC -**
075014	Power supply charger 220 V	**NC -**
075016	Complete culture nosepiece with 4 filters	**NC -**
075131	Complete valve + 2 buttons for Pipet-Aid XP	**NC -**
075132	Complete valve + 2 buttons for Pipet-Aid XP Gravity	**NC -**
075133	Valve only for Pipet-Aid	**NC -**
075776	Hydrophobic filter 0,2 µm (5 units)	**NC -**
	Spare parts for Pipet-Aid XP2	
075214	Filter lock	**NC -**
075215	Filter connections	**NC -**
075216	Removable stand for Pipet-Aid XP2	**NC -**
075195	Batteries	**NC -**

Pipet-Aid XP2 : Even better performance!



