

EasyFit® Treff® Microtubes

	CleanRoom Pure®					LoTrace® DNA	
EasyFit®							
Volume (ml)	0.5		1.5		2.5	0.5	1.5
Diameter (mm)	10	10	12.8	12.8	12.8	10	12.8
Height (mm)	30	30	38.8	38.8	38.7	30	38.8
Material	Virgin polypropylene	Virgin polypropylene	Virgin polypropylene	Virgin polypropylene	Virgin polypropylene	Virgin polypropylene	Virgin polypropylene
Sterile	No	Yes	No	Yes	No	No	No
Cap Type	Snap EasyFit®						
Autoclavable	Yes						
Max centrifugal speed (g)	30000						
Graduation	Yes						
Marking area	Yes						
Operating temperature	-195 °C to +121 °C						
Specifications	For general use; easy opening; cleanroom class 8 according to ISO 14644-1; free of DNases, RNases, human DNA, endotoxins, ATP and PCR inhibitors					Low DNA adhesion; for general use; easy to open; cleanroom class 8 according to ISO 14644-1; free of DNases, RNases, human DNA, endotoxins, ATP and PCR inhibitors	
Units/box	500 microtubes		1000 microtubes			2000 microtubes	1000 microtubes
Natural	**034001**	**034159**	**034013**	**034013S**	**982854**	**982748**	**982749**
Yellow	**034003**	-	**034015**	-	-	-	-
Green	**034004**	-	**034016B**	-	-	-	-
Blue	**034005**	-	**034017**	-	-	-	-
Pink	**034002**	-	**034014**	-	-	-	-
€ Excl. VAT/box	**NC -**	**NC -**	**NC -**	**NC -**	**NC -**	**NC -**	**NC -**

TIP





➤CapLock :
Fully sealed caps, double lock to prevent opening of cap during incubation



➤EasyFit :
For general use and storage



➤ClickFit :
Reinforced closure, particularly suitable for heating applications