Bouchons gris en mousse ClearLine® Grey foam caps



- In polyether foam P2560
- Reusable
- Smooth surface for marking with a felt-tip pen
- Density 25 kg/m³
- Support 6.3 kPa

Cat. No.	Diameter (mm)	Height (mm)	Units/ carton	€ Excl. VAT/ box
999046CL	50	50	100	NC -
999047CL	40	40	100	NC -
999048CL	26	30	100	NC -

Cotton or rayon balls and ropes



- For capping Drosophila flasks and tubes
- The ropes are cut to the selected dimensions

Cat. No.	Description	Units/ carton	€ Excl. VAT box
789041	Rayon ball for 25 mm tube	2000	NC -
789042	Cotton ball for 25 mm tube	2000	NC -
789043	Cotton ball for 28.5 mm tube	2000	NC -
789044	Cotton ball for bottle	2000	NC -
789301	Cotton ball for extra large bottle	1000	NC -
789302	Rayon ball for extra large flask	2000	NC -
789045	Roll of rayon wire	7.5 kg	NC -
789046	Roll of cotton wire	7,5 kg	NC -

White or grey foam caps



- Reusable
- White or Grey
- Smooth surface for marking with a felt-tip pen

Cat. No.	Description	Colour	Units/ carton	€ Excl. VAT box
330070	Foam cap ø x H = 36 x 40 mm	White	605	NC -
999036	Foam cap ø x H = 50 x 50 mm	Grey	100	NC -

Gaze Cheese cloth



- Made of 100% woven cotton
- 28 x 24 wires per square inch
- Unfolded : 91.4 cm x 64 m

Cat. No.	Description	€ Excl. VAT
789056	Gaze Cheese cloth	NC -

Filter mesh

Made of stainless steel or nylon.



Stainless steel filter mesh

Cat. No.	Description	€ Excl. VAT
789048	Stainless steel filter mesh 97 µm	NC -

NITEX nylon filter cloths

Cat. No.	Porosity (µm)	Applications	€ Excl. VAT
074002	500	Adult retrieval	NC -
074004	250	Various	NC -
074007	100	Embryo retrieval	NC -
074010	50	CO_2 Filters and nuclear debris filters during the preparation of the larvae	NC -
074013	30	Various	NC -

GeneSeal adhesive film



Cat. No.	Description	€ Excl. VAT
789163	GeneSeal film for 24 mm ø - roll 192 m	NC -
789164	GeneSeal film for 28.5 mm ø - roll 192 m	NC -

Droso-Filler filling device



Food filling device for tubes.

- To fill 100 thin diameter tubes
- To fill 100 large diameter tubes
- To fill 25 square bottles at once using the tray

Cat. No.	For	Suitable tubes (mm)	Ext. dimensions (cm)	Int. dimensions (cm)	Distance between holes (cm)	€ Excl. VAT/ box
789166	Thin tubes	Ø 25	31.4 x 34.9 x 29.2	26 x 27.9 x 14	2.54	NC -
789167	Large tubes	Ø 28,5	35.5 x 36.8 x 30.5	29.8 x 30.2 x 14	3.17	NC -
789168	Square bottles	57 x 57	35.5 x 36.8 x 30.5	29.8 x 30.2 x 14	5.71	NC -