PCR refrigerated racks

- Polypropylene racks, polycarbonate lid
- Dimensions (cm) : 11 x 14.8 x 4.8
- Compatible with PCR tube strips or individual 0.2/0.5 ml PCR tubes, in standard or low profile format
- Maintains the temperature of a 96-well PCR plate at 0°C (after cooling the rack to -20°C for 2 hours) for more than 1 hour
- Maintains the temperature of the samples below +5°C for up to 3 hours
- Change of colour of the rack or a well when the temperature exceeds + 6 °C
- The lid provides optimum protection for samples from the freezer to the bench
- Innovative, sliding, 180° rotating lid (can also be used as a stand)
- Suitable for storage at -80 °C
- Wells designed to prevent freezing of samples



320688



320689

Cat. No.	Colour non-refrigerated	Colour refrigerated	Units/ box	€ Excl. VAT/box
320688	Pink	Violet	2	NC -
320689	Light blue	Blue	2	NC -
320690	1 rack of each colour		2	NC -

Eppendorf IsoTherm-System

eppendorf

- Insulated box and rack system with cold storage blocks to maintain samples at 0°C or -21°C during bench top handling
- Ideal for enzymatic workstations, PCR stations or sample transfer
- Replaces the ice bath
- Compatible with the accessories of the Thermomixer range
- The system includes:
  - an IsoSafe box accepting an IsoRack and an IsoPack cold accumulator
  - an IsoRack for 24 microtubes, microplate format, centrifugal, floating
  - an IsoPack cold accumulator: 0 °C with 6h runtime (white), or -21 °C with 3h runtime (blue)



For microtubes	1.5 / 2 ml				0.5 ml			
Description	Cat. No. 0 °C	€ Excl. VAT	Cat. No. -21 °C	€ Excl. VAT	Cat. No. 0 °C	€ Excl. VAT	Cat. No. -21 °C	€ Excl. VAT
1 IsoRack + 1 IsoPack	033579	NC -	033612	NC -	033613	NC -	033614	NC -
1 IsoSafe + 3 IsoPack	033615	NC -	-	-	033617	NC -	033618	NC -
	Cat. No.		€ Excl. VAT		Cat. No.		€ Excl. VAT	
IsoRack with lid	033619		NC - 4 units		033621		NC - 4 units	
Complete system : 1 Isosafe + 1 IsoPack 0 °C + 1 IsoPack -21 °C +1 IsoRack with lid	033620		NC -		033622		NC -	

CoolCADDY cold maintenance Workstation

- Temperature + 2 °C - + 4 °C up to 5 hours
- Keeping samples cold and unfrozen
- Suitable for :
  - 9 x 1.5/2 ml microtubes;
  - 8 x 2 ml cryotubes ;
  - 5 x 0.5 ml microtubes ;
  - 1 plate;
  - 96 PCR microtubes or 12 strips;
  - 1 x 50 ml tank



Cat. No.	Description	€ Excl. VAT
062900	CoolCADDY cold maintenance Workstation	NC -

Reagent cooling tanks

- Made of polypropylene
- Ice-free, non-toxic
- To increase handling time while preserving the integrity of sensitive reagents
- Reagent held at less than +5°C for approximately 3 hours
- Suitable for tanks Cat. No. 320661 to 320669
- Stackable
- Reusable (warranty : 2 years)



Vol. (ml)	Dimensions (mm)	Pink	Blue	Assorted	Units/ box	€ Excl. VAT/ carton
50	144 x 74 x 45	320695	320696	320697	2	NC -
100	150 x 90 x 43	320698	320699	320700	2	NC -

Indicative prices. Visit www.dulis.nl