Single-use Scissors

- **Stainless steel blade and plastic handle**
- **Sterilised with ethylene oxide**
- **Unit packaging**
- **Pack with 100 units**

Cat. No.	Length (mm)	Ends	Type	Units/ box	€/ box
030676	115	Rounded	-	100	**NC -**
030677	115	Pointed/Foam	High accuracy	100	**NC -**
030678	115	Pointed	-	100	**NC -**
030679	110	Pointed	Wire Cutter	100	**NC -**


030678


030677


030676


030679

Wironit 18/12 stainless steel scissors

Wironit is an 18/12 stainless steel with 18% chromium and 10-12% nickel, which is even more corrosion resistant and anti-magnetic.


491010


491011


491012


491013


491014

Cat. No.	Length (mm)	Extremities	€
491010	110	Pointed - Pointed	**NC -**
491011	145	Pointed - Pointed	**NC -**
491012	130	Pointed - Rounded	**NC -**
491013	145	Pointed - Rounded	**NC -**
491014	160	Pointed - Rounded	**NC -**

Universal scissors



All stainless steel universal scissors

Cat. No.	Type	Total Length (mm)	Blade length (mm)	€
076139	All stainless steel	100	30	**NC -**
076140	All stainless steel	115	45	**NC -**
076141	All stainless steel	130	50	**NC -**
076142	All stainless steel	145	60	**NC -**
076144	All stainless steel	165	70	**NC -**
076145	All stainless steel	180	80	**NC -**



All-purpose scissors in stainless steel, particularly suitable for the laboratory.

Universal scissors with plastic handle

Cat. No.	Type	Total Length (mm)	Blade length (mm)	€
076146	Plastic handle	130	50	**NC -**
076147	Plastic handle	150	70	**NC -**
076148	Plastic handle	200	100	**NC -**
076149	Plastic handle	250	120	**NC -**

All prices are indicative and exclude VAT. Visit www.dulis.nl