
Glass Centrifuge Tubes, graduated

- Cylindro-conical
- In AR-GLAS or DURAN® glass
- Capacity: approx. 15 ml
- Inscriptions and graduations in high contrast white enamel



Cat. No.	Grad. (ml)	Glass	ø x H (mm)	Subdivisions (ml)	Thickness (mm)	ACR max.	€ 10 units
092305	0 - 10	AR-GLAS	17 x 113	0.1	1	3500	NC -

Centrifuge tubes, thick wall, screw cap

- Round bottom
- With screw cap with PTFE rubber seal
- Borosilicate glass
- Thick wall for resistance to higher speeds
- Autoclavable tubes and screw caps



Cat. No.	Volume (ml)	ø x H (mm)	Units/ box	€/ box
045608	35	29 x 100	12	NC -
045609	50	29 x 123	12	NC -

Replacement cap

Cat.No.	Description	Units/ box	€/ box
045610	Replacement caps	150	NC -

High-Speed Glass Tubes, super-resistant

- Reusable
- Made of special borosilicate glass type I, class B
- Chemically reinforced for improved mechanical strength
- Can be centrifuged up to 13100 RCF in a 50 ml tube cavity with the adapter
- Can withstand temperatures up to 300°C
- With straight neck without cap, or with screw cap with PTFE filmed rubber seal

Volume (ml)	15	30
ø x H (mm)	18 x 102	24 x 106
Tubes with straight neck		
Cat. No.	045600	045601
€/ 6 units	NC -	NC -
Screw-cap tubes		
Cat. No.	045602	045603
€/ 6 units	NC -	NC -
Rubber adapters for 50 ml rotor		
Cat. No.	045604	045605
€/ 2 units	NC -	NC -
Replacement caps for screw-in tubes		
Cat. No.	-	045243
€/ 144 units	-	NC -



Tubes for Hybridization

Borosilicate glass tubes for Northern, Southern or Western blots hybridization. With watertight polypropylene screw caps with silicone O-Ring diameter 45 mm.

Cat. No.	ø x W (mm)	Units/ box	€/ box
045215	35 x 150	2	NC -
045216	35 x 300	2	NC -
001134	40 x 170	1	NC -
001135	40 x 260	1	NC -
001136	75 x 260	1	NC -



Find more hybridization tubes on dulis.nl



All prices are indicative and exclude VAT. Visit www.dulis.nl