Pig HAZ-MAT Chemical Absorbent pads


333202


333206



333207


333205


333208

- 8 layers of fine polypropylene fibres
- For the absorption of a wide range of chemicals (e.g. 98% sulphuric acid, 30% sodium hydroxide)
- Pads with honeycomb pattern: faster absoption

Pads in a bag

Cat. No.	Description	Dimensions (mm)	Absorption (l/bag)	Pads/ bag	€ the bag
333202	Heavyweight mat pads	380 x 510	41.8	50	NC -
333203	Heavyweight mat pads	380 x 510	83.6	100	NC -
333204	Lightweight mat pads	380 x 510	83.6	200	NC -

Pads in a box

Cat. No.	Description	Dimensions (mm)	Absorption (l/box)	Pads/ box	Box/ cardboard box	€/ cardboard box
333206	Heavyweight mat pads in Box dispenser,	250 x 330	37.1	100	1	NC -
333207	Heavyweight mat pads in Bench Box dispenser,	330 x 250	37.1	25	4	NC -

Roll

Rip&fit rolls: delivered in a dispenser box.

Cat. No.	Description	Dimensions (mm x m)	Absorption (l/bag)	Rolls/ bag	€ the bag
333205	Heavyweight mat rolll	760 x 46	152.4	1	NC -
333208	Heavyweight mat roll in Rip&fit dispenser box	380 x 18	29.7	1	NC -

Lab Sorb

- For gelling of liquid products
- Lab Sorb granules gel aqueous compounds or solvents in minutes
- Gelled product is easily disposed of
- 450 g



Cat. No.	Description	€
053489	Lab Sorb, 450 g	NC -


Universal absorbent for oils and chemicals



- Absorbs chemicals or oils
- Fine polymer-based granules, dust-free
- For all liquids
- Absorbs up to 75 times its own weight
- Economic
- Coloured indicator depending on the compound to be absorbed: yellow for acid-based products, red for alkaline products
- Simply sprinkle on the spilled product, collect in gel form with a shovel or spatula (not supplied)

Cat. No.	Description	€
452000	Chemsolute chemical absorbent, 1800 ml	NC -
452001	Chemsolute chemical absorbent, 4000 ml	NC -

HandySorb Absorbents



- Inert polypropylene sheets and cushions
- Each leaf absorbs up to 25 times its weight
- Do not react with any spilled substances
- NoTouch system: no user contact with the spilled liquid
- User range: up to 1.83m for more effective use in hard-to-reach areas

Cat. No.	Description	Units/ box	€/ box
834423	HandySorb System Starter Kit with absorbent pads	1	NC -
834420	Pads refills with adapter for HandySorb system	100	NC -
834421	NoTouch universal pads for HandySorb system	25	NC -
834422	NoTouch universal cushions for HandySorb system	10	NC -

Beta-safe paper



- Activated carbon paper to absorb volatile radioisotopes (^{35}S)
- Prevents contamination of laboratories

Cat. No.	Description	Units/ box	€/ box
034918	Beta-Safe sheets 20.5 x 25.4 cm	50	NC -