Pipetting device





Pipetting device for all volumetric and graduated pipettes from 0.1 to 200 ml.

- Suction: once the bulb is compressed, the lever is operated to suck up the liquid, with easy speed control.
- Distribution: the liquid is discharged by lowering the lever.
- Purge: By pressing the rubber button, the last drops (meniscus) are ejected from the pipette.

Blue	Pink	Green	Grey	€ Excl. VAT
092475	092476	092474	092473	NC -

Accessories

Cat.No	Description	Units/ box	€ Excl. VAT/ box
092480	Filter membrane 3 µm (PP, PTFE), non-sterilised	1	NC -
042350	Filter membrane 3 µm (PP, PTFE), non-sterilised	10	NC -
092379	Adapter (silicone), 44 mm length	1	NC -
092478	Adapter box (PP), 49 mm length	1	NC -
092479	Suction bellows (silicone), screw ring (PP)	1	NC -
092481	Valve system (PP, PTFE, silicone)	1	NC -

Blue pipetting device



- **Dimensions (mm) 230 x 54 x 54**
- **Ergonomic handle**
- **Compatible with plastic or glass serological pipettes from 1 to 100 ml**
- **Control of the suction/distribution speed by means of a T-bar**
- **2 extra tips (green and blue) included**
- **45 µm hydrophobic filter**
- **PP tip and silicone adapter are autoclavable**

Cat.No.	Description	€ Excl. VAT
320554	Blue pipetting device	NC -

Pipetting aid



- **Comfortable and easy to use**
- **Precise pipetting control**
- **Robust and lightweight**
- **Easy maintenance**
- **Suitable for volumetric and graduated pipettes, glass or plastic, from 0.1 to 100 ml**
- **Replaceable 3 µm Spare filter**

Colour	Blue	Red	Green	Yellow	€ Excl. VAT
Cat. No.	861000B	861001B	861002B	861003B	NC -

Accessory

Cat. No.	Description	€ Excl. VAT
861005	Spare filter 3 µm	NC -