POLYCARBONATE/POLYPROPYLENE TUBES, STRAIGHT NECK

Nalgene™ Polycarbonate Tubes, straight neck

Thermo SCIENTIFIC NALGENE



Volume (ml)	12	15		16	38		50		100
out.Ø (mm)	15.9	17	16	17.9	25.3	28.7	28.8	28.5	31.5
Height (mm)	103.1	120.1	114	99.8	88.6	134.1	103.1	102	164.3
Material	Polycarbonate								
Sterile	No								
Bottom	Rounded	Conical	Rounded	Rounded	Rounded	Conical	Rounded	Rounded	Rounded
Caps	Without								
Type of neck	Straight	Straight	Straight	Straight	Straight	Straight	Straight	Straight with edge	Straight with edge
Autoclavable	Yes								
Max centrifugal speed (g)	50000	6000	50000	50000	50000	6000	50000	50000	50000
Skirt	No								
Graduation	No	Yes	No	No	No	Yes	No	No	No
Marking area	No								
Colour	Transparent								
Operating temperature	-135 °C to +135 °C								
Units/box	Bag of 10 tubes								
Cat. No. tube only	**029222**	**029209**	**029223**	**029224**	**029225**	**029210**	**029226**	**029228**	**029227**
€ Excl. VAT/ box	**NC -**	**NC -**	**NC -**	**NC -**	**NC -**	**NC -**	**NC -**	**NC -**	**NC -**

Polypropylene Tubes, straight neck

Round bottom
Without cap



Volume (ml)	1.6	3	3.1	5								
out.Ø (mm)	8	11.5	12	12	12.8	13	12	12	13	13	13	12
inn.Ø (mm)	-	-	-	-	10.9	11	-	-	-	-	-	-
Height without cap (mm)	50	55	50	75	75	75	75	75	75	75	75	86
Material	Polypropylene											
Sterile	No											
Bottom	Rounded											
Type of neck	Straight	Straight	Straight	Straight	Straight	Straight	Straight	Straight	Straight	Straight	Straight with edge	Straight with edge
Autoclavable	-	No	-	-	No	No	No	No	No	No	No	No
Max centrifugal speed (g)	10000	5000	10000	10000	-	-	5000	5000	5000	5000	5000	5000
Skirt	No											
Graduation	No	No	No	No	No	No	No	No	No	No	Yes	Yes
Marking area	No											
Colour	Translucent											
Specifications	-	-	-	-	-	-	Sorvall tube	Sorvall tube	-	-	-	-
Units/box	Bag of 1000 tubes	12 bags of 500 tubes	Bag of 1000 tubes	Bag of 1000 tubes	Bag of 2000 tubes	4 bags of 1000 tubes	16 bags of 250 tubes	Bag of 4000 tubes	16 bags of 250 tubes	Bag of 4000 tubes	Bag of 4500 tubes	Bag of 4000 tubes
Cat. No. tube only	**064375**	**390017**	**064317**	**064200**	**080071**	**388589**	**390018**	**390019**	**390020**	**390021**	**390026**	**390027**
€ Excl. VAT/ box	**NC -**	**NC -**	**NC -**	**NC -**	**NC -**	**NC -**	**NC -**	**NC -**	**NC -**	**NC -**	**NC -**	**NC -**

Indicative prices. Visit www.dulis.nl