CAPS

Snap caps with tab

- Made of polyethylene
- For plastic or glass tubes
- Easy to open thanks to the grip tab



ø (mm)	Tab	Grey	Blue	Natural	Green	Lavender	Red	White	Yellow	€ Excl. VAT 1000 uni.
12 - 13	1	380527	380528	380529	380530	380531	380532	380533	380534	NC -
16	1	380535	380536	380537	380538	380539	380540	380541	380542	NC -
12	2	-	380596	380597	380598	380599	380600	-	380601	NC -
13	2	380602	380603	380604	380605	380606	380607	-	380608	NC -
16	2	380609	380610	380611	380612	380613	380614	-	380615	NC -

Drillable caps

- The material used does not interfere with common chemicals or coagulation
- Can be used at room temperature, at 4°C or in the freezer at -20°C
- For tube with outer Ø 12-13 mm or 16 mm
- In santoprene
- Single use
- Can be applied with one hand
- Can be drilled with a needle



Colour	For 12 - 13 mm tubes	For 16 mm tubes
Blue	040245	040252
Green	040246	040253
Gray	040247	040254
Lavender	040248	040255
Red	040249	040256
White	040250	040257
Yellow	040251	040258
€ Excl. VAT 1000 units	NC -	NC -

Autoclavable vented caps

Polypropylene caps, autoclavable, with internal ventilation lug. Reusable



Tube ø (mm)	Units/ box	Natural	€ Excl. VAT box	Red	Green	Yellow	Blue	Units/ box	€ Excl. VAT box
13	500	045082	NC -	045088	045094	065100	065105	500	NC -
16	500	045083	NC -	045089	045095	065101	065106	500	NC -
16	1000	140434	NC -	140436	140439	140437	140438	1000	NC -
18	500	045084	NC -	045090	045096	065102	065107	500	NC -
20	250	045085	NC -	045091	045097	065103	065108	250	NC -
25	250	045086	NC -	045092	045098	065104	065109	250	NC -
38	125	045087	NC -	-	-	-	-	-	-

Double position caps

Simport

- For bacteriology tubes
- Closed position: anaerobic environment / open position: aerobic environment
- Made of polyethylene

Cat. No.	Double position cap	Units/ box	€ Excl. VAT box
039387	For polystyrene tube 12 mm Ø	1000	NC -
039388	For polypropylene tube 12 mm Ø	1000	NC -
039269	For polystyrene tube 17 mm Ø	1000	NC -
039389	For polypropylene tube 17 mm Ø	1000	NC -

Indicative prices. Visit www.dulis.nl