
FLOATING RACKS/ FOR BOTTLES

Floating foam racks



Foam racks for tubes and microtubes floating in a water bath, with handle.

Cat. No.	Description	Colour	€ Excl. VAT 5 units
037544	Racks for 4 x 50 ml tubes, 30 mm ø	Green	NC -
037545	Racks for 8 x 15 ml tubes, 17 mm ø	Blue	NC -
037546	Racks for 24 x 1,5 ml or 2 ml microtubes	Yellow	NC -

Floating rack for PCR microtubes



- In red foam
- For 10 PCR tubes or for 4 x 8-microtube strips

Cat. No.	Description	€ Excl. VAT 10 units
030122	Floating rack for 10 PCR tubes	NC -
030123	Floating rack for 4 x 8-PCR tube strips	NC -

Floating rack for centrifuge tubes



- Protects samples from loss and contamination
- For 28 x 1.5 to 2 ml centrifuge tubes
- Autoclavable jup to 121 °C

Cat. No.	Description	€ Excl. VAT
086700	Floating rack for centrifuge tubes	NC -

Floating foam rack for microtubes



Floating foam rack to fit in 500 ml beakers.

- From 18 microtubes from 0.2 to 1.5 - 2 ml
- With handle to be fixed
- Diameter : 7.62 cm

Cat. No.	Description	€ Excl. VAT 5 units
037594	Floating foam rack for 18 microtubes	NC -

Floating racks for microtubes



- Made of polypropylene
- Tab for holding the tube to the rack
- Carrying handles
- Autoclavable

Cat. No.	Capacity	Dimensions W x D x H (mm)	Units/ box	€ Excl. VAT/box
947717	60 x 1.5- 2 ml microtubes	269 x 77 x 43	2	NC -
088899	24 x 5 ml microtubes	170 x 79 x 54	2	NC -

Holder for bottles or round flasks



- For different volumes of flasks
- Excellent resistance to organic solvents
- Operating temperature range : - 30 to + 120 °C

Cat. No.	Number of flasks	Volume (ml)	Hole ø x H (mm)	W x D x H (mm)	€ Excl. VAT
446711	3	100	70 x 91	310 x 100 x 100	NC -
446706	12	100	70 x 91	384 x 265 x 98	NC -
446712	2	250	98 x 122	280 x 125 x 130	NC -
446707	6	250	98 x 122	389 x 239 x 130	NC -
446708	2	500	121 x 128	360 x 146 x 135	NC -
446709	2	1000	148 x 156	415 x 192 x 165	NC -
446710	2	2000	183 x 188	503 x 223 x 199	NC -

Collection baskets



- Side straps
- Dimensions : 432 x 279 x 279 mm

Colour	Blue	€ Excl. VAT
Cat. No.	320119	NC -