High density polyethylene bottles, wide neck

Nalgene™ round bottles with wide neck

Thermo SCIENTIFIC NALGENE

- Round bottle with wide neck in HDPE
- Waterproof
- Polypropylene cap





Capacity (ml)	30	60	125	250	500	1000	1500
Cap size (mm)	28	28	38	43	53	63	63
Units/bag	12	12	12	12	12	6	6
Units/box	72	72	72	72	48	24	6
HDPE Round bottle with wide neck, PP cap							
Cat. No.	**028725**	**028726**	**028727**	**028728**	**028729**	**028730**	**029391**
€/box	**NC -**	**NC -**	**NC -**	**NC -**	**NC -**	**NC -**	**NC -**
Amber HDPE round bottle with wide neck, PP amber cap							
Cat. No.	**028744**	**028745**	**028746**	**028747**	**028748**	**028749**	**-**
€/box	**NC -**	**NC -**	**NC -**	**NC -**	**NC -**	**NC -**	**-**

Nalgene™ economy bottles



HDPE bottle with wide neck, PP cap

Thermo SCIENTIFIC NALGENE

Features	Waterproof					
Capacity	30	60	125	250	500	1000
Cap size (mm)	28	28	38	43	53	63
Cat. No.	**028810**	**028811**	**028812**	**028813**	**028814**	**028815**
Units/box	72	72	72	72	48	24
€/box	**NC -**	**NC -**	**NC -**	**NC -**	**NC -**	**NC -**
Cat. No.	**-**	**-**	**028820**	**028821**	**028822**	**028823**
Units/box	-	-	500	250	125	50
€/box	**-**	**-**	**NC -**	**NC -**	**NC -**	**NC -**
Amber HDPE bottle with wide neck, amber PP cap						
Features	Waterproof					
Cat. No.	**-**	**-**	**028806**	**028807**	**028808**	**028809**
Units/box	-	-	72	72	48	24
€/box	**-**	**-**	**NC -**	**NC -**	**NC -**	**NC -**

HDPE bottles

- In translucent HDPE
- Polypropylene screw caps
- Waterproof
- Good chemical resistance

DWK LIFE SCIENCES



		10 units		BULK version		
Vol. (ml)	**Ø x H (mm)**	**Cat. No.**	**€**	**Cat. No.**	**Units/ box**	**€/ box**
30	32 x 70	**723861**	**NC -**	**761838**	1200	**NC -**
60	40 x 82	**723862**	**NC -**	**761839**	740	**NC -**
100	51 x 89	**232659**	**NC -**	**761840**	420	**NC -**
125	51 x 104	**232581***	**NC -**	**761841***	350	**NC -**
150	51 x 116	**232660***	**NC -**	**761842***	350	**NC -**
250	60 x 140	**723883**	**NC -**	**761843**	200	**NC -**
500	75 x 161	**723864**	**NC -**	**761844**	100	**NC -**
1000	90 x 205	**723865****	**NC -**	**761845**	60	**NC -**
2000	120 x 242	**121304****	**NC -**	**761846**	30	**NC -**

* Not approved for food contact ** Per 5

Wheaton® HDPE bottles with wide neck

DWK LIFE SCIENCES

- Waterproof
- Anti-drip ring
- Delivered with cap assembled
- Double-sealed cap



Volume (ml)	30	60	125	250	500	1000
ø x H (mm)	33 x 69	39 x 88	51 x 98	61 x 125	72 x 164	91 x 199
Cap	28-410	28-410	38-410	43-410	53-410	63-415
Units/bag	12	12	12	12	12	12
Units/box	72	72	72	72	48	48
HDPE bottle						
Cat. No.	**048353**	**048354**	**048355**	**048356**	**048357**	**048358**
€/bag	**NC -**	**NC -**	**NC -**	**NC -**	**NC -**	**NC -**
Amber HDPE bottle						
Cat. No.	**048369**	**048360**	**048361**	**048362**	**048363**	**048772**
€/bag	**NC -**	**NC -**	**NC -**	**NC -**	**NC -**	**NC -**

All prices are indicative and exclude VAT. Visit www.dulis.nl