Laboplast and Steriplast sampling scoops, spoons and spatulas







- White polystyrene
- Individual packaging
- Produced and assembled in a class 7 clean room
- Steriplast quality: sterile gamma irradiated (minimum 10 kGy)
- Complies with EU food and FDA directives

Laboplast bio and Steriplast bio

- White polyethylene
- Materials and packaging made entirely from renewable resources

Steriplast detectable

- In blue polystyrene
- Detectable by colour scanner, metal detector or X-ray

GoGREEN

1 Sampling scoops

Capacity (ml)	Length (mm)	Lid	Units/ bag	Units/ box	Laboplast Cat. No.	Laboplast €	Steriplast Cat. No.	Steriplast €	Laboplast bio Cat. No.	Laboplast bio €	Steriplast bio Cat. No.	Steriplast bio €	Steriplast detectable Cat. No.	Steriplast detectable €
25	141	No	1	10	795500	NC -	795508	NC -	795516	NC -	795520	NC -	795651	NC -
25	142	Yes	1	10	795501	NC -	795509	NC -	-	-	-	-	-	-
50	160	No	1	10	795502	NC -	795510	NC -	795517	NC -	795521	NC -	795652	NC -
50	161	Yes	1	10	795503	NC -	795511	NC -	-	-	-	-	-	-
100	205	No	1	10	795504	NC -	795512	NC -	795518	NC -	795522	NC -	795653	NC -
100	206	Yes	1	10	795505	NC -	795513	NC -	-	-	-	-	-	-
150	216	No	1	10	795506	NC -	795514	NC -	795519	NC -	795523	NC -	795654	NC -
150	218	Yes	1	10	795507	NC -	795515	NC -	-	-	-	-	-	-
250	232	No	1	10	795560	NC -	795564	NC -	-	-	-	-	795655	NC -
250	235	Yes	1	10	795561	NC -	795565	NC -	-	-	-	-	-	-
500	280	No	1	10	795562	NC -	795566	NC -	-	-	-	-	795656	NC -
500	283	Yes	1	10	795563	NC -	795567	NC -	-	-	-	-	-	-
1000	332	No	1	10	795639	NC -	795640	NC -	-	-	-	-	795657	NC -

2 Sampling spoons

Capacity (ml)	Length (mm)	Units/ bag	Units/ box	Laboplast Cat. No.	Laboplast €	Steriplast Cat. No.	Steriplast €	Laboplast bio Cat. No.	Laboplast bio €	Steriplast bio Cat. No.	Steriplast bio €	Steriplast detectable Cat. No.	Steriplast detectable €
2.5	127	1	100	795524	NC -	795526	NC -	795528	NC -	795530	NC -	795659	NC -
10	170	1	100	795525	NC -	795527	NC -	795529	NC -	795531	NC -	795660	NC -

3 Sampling spatulas

L. blade (mm)	L. total (mm)	Sheath	Units/ bag	Units/ box	Laboplast Cat. No.	Laboplast €	Steriplast Cat. No.	Steriplast €	Laboplast bio Cat. No.	Laboplast bio €	Steriplast bio Cat. No.	Steriplast bio €	Steriplast detectable Cat. No.	Steriplast detectable €
150	263	No	1	10	795532	NC -	795534	NC -	795536	NC -	795537	NC -	795658	NC -
150	268	Yes	1	10	795533	NC -	795535	NC -	-	-	-	-	-	-





SteriPlast Sterile Sampling Kit

Kit containing ten Steriplast scoops or spatulas and ten SteriBag Premium bags.

Cat. No.	Kit components	€
795662	10 Steriplast 50 ml sampling scoops Cat. No. 795510 and 10 SteriBag Premium 300 ml Cat. No. 795661	NC -
795663	10 Steriplast 150 mm sampling spatulas Cat. No. 795532 and 10 SteriBag Premium 300 ml ref. 795661	NC -

Small spatulas



- White
- Withstands down to 230°C

Cat. No.	Length (cm)	€
857049	30	NC -
857050	40	NC -



- 200 °C resistant handle
- Head at 260 °C

Cat. No.	Length (mm)	€
857051	250	NC -
857052	350	NC -
857053	450	NC -