
Nunc CryoBank™ and Bank-It™ Storage tubes

For automated storage of samples. Volume 0.5 ml or 1 ml

- **CryoBank™ for cells and proteins**
- **Bank-It™ for DNA**
- **Optimised storage: microplate format, storage in a cryobox**
- **96 microtubes in 8 by 12 format on microplate, or in bulk**
- **The tubes are clicked into the rack: easy opening and closing**
- **Low binding surface for improved retrieval of low concentration samples**
- **100% leak-proof double click cap**
- **2D barcode, 14 x 14**



	2D Barcode	CryoBank 0,5 m Tubel	CryoBank 1 ml Tube	Bank-It 0,5 m Tubel	Bank-It 1 ml Tube
Tubes on rack	No	**056282**	-	-	-
€ Excl.VAT the 2 x 96		**NC -**	-	-	-
Tubes on rack	Yes	-	**056285**	-	-
€ Excl.VAT the 2 x 96		-	**NC -**	-	-
Tubes on rack	Yes	**056289**	**056293**	**056291**	**056295**
€ Excl. VAT the 10 x 96		**NC -**	**NC -**	**NC -**	**NC -**
Tubes Sl in bulk	No	**056297**	**056298**	-	-
€ Excl. VAT the 10 x 96		**NC -**	**NC -**	-	-

Accessories

Cat. No.	Description	Units/ box	€ Excl. VAT/ box
056364	Freezing shelf (H=74 cm) for 19 boxes of 0.5 ml tubes + 13 boxes of 1 ml tubes	1	**NC -**
056365	Single caps for Bank-It cryotubes	960	**NC -**
056540	Green caps for Bank-It cryotubes	960	**NC -**
056541	Red caps for Bank-It cryotubes	960	**NC -**
056542	Blue caps for Bank-It cryotubes	960	**NC -**
056299	Electric screwdriver	1	**NC -**
056300	Manual screwdriver	1	**NC -**

Nunc Cryobank 2 and 5 ml Cryotubes


- **2 or 5 ml screw cap tubes on sterile rack of 48 tubes or in bulk, with or without 2D code**
- **Rack format to SBS standard**
- **Self-standing tubes, self-locking in the rack**
- **Same dimensions as conventional cryotubes**
- **Translucent**



Tot. Vol. µl/wells	2 ml	5 ml	2 ml	5 ml
2D Code	With	With	Without	Without
On rack	**056733**	**056732B**	**056735**	**056734**
Units/carton	48 x 10	48 x 10	48 x 10	48 x 10
€ Excl. VAT/carton	**NC -**	**NC -**	**NC -**	**NC -**
Bulk	**056775**	**056779**	**056777**	**056778**
Units/carton	960	528	960	528
€ Excl. VAT/carton	**NC -**	**NC -**	**NC -**	**NC -**
On lockable rack	**056797**	**056796**	-	-
Units/carton	48 x 10	48 x 10	-	-
€ Excl. VAT/carton	**NC -**	**NC -**	-	-

Cat. No.	Description	€ Excl. VAT/ 10 units
056792	Rack for 2 ml Cryobank tubes	**NC -**
056793	Rack for 5 ml Cryobank tubes	**NC -**

Nunc "Press out tool"



Tool for rapid retrieval of Bank-It and CryoBank cryotubes Allows you to select the tubes to be removed from a complete rack.
Use :

- place the metal pegs on the selected position on the tool;
- insert the rack with the tubes on the tool (the orientation is controlled by the cut-off edge of the rack);
- with both handles, press the rack, the selected tubes are removed from the rack. Put the rack back in the cold. The operation took only a few seconds and the samples do not have time to thaw.

Cat. No.	Description	€ Excl. VAT
056330	Nunc "Press out tool"	**NC -**

2D transfer tubes 1.8 ml, 2 ml and 5 ml

- **Polypropylene tube, polycarbonate rack**
- **USP Class VI**
- **Cap with external polypropylene thread with thermoplastic elastomer seal**
- **Low version, total tube height 49.79 mm**
- **Temperature resistant from -180 °C to 121 °C**
- **Sterile, autoclavable**
- **Graduated, white marking area**



	5 ml tubes				2 ml tubes with internal screw thread		1.8 ml tubes	
Type/ Marking	**Cat.No**	**Units/ carton**	**€ Excl. VAT/ carton**	**Units/ carton**	**Cat.No**	**€ Excl. VAT/ carton**	**Cat.No**	**€ Excl. VAT/ carton**
Latch rack/2D	**056851**	24 x 10	**NC -**	48 x 10	**056872**	**NC -**	**056843**	**NC -**
Latch rack	**056852**	24 x 10	**NC -**	48 x 10	**056873**	**NC -**	**446861**	**NC -**
In bulk/2D	**056853**	240	**NC -**	480	**056874**	**NC -**	**056844**	**NC -**
In bulk	**056854**	240	**NC -**	480	**056875**	**NC -**	-	-
Empty rack	**056855**	10	**NC -**	10	**056871**	**NC -**	**262181**	**NC -**
Empty latch rack with barcode	-	-	-	10	-	-	**056906**	**NC -**
Latch rack with barcode	-	-	-	48 x 10	**056872B**	**NC -**	**056843B**	**NC -**

Indicative prices. Visit www.dulis.nl