
Nunc™ plates


Omnitray Nunc™ culture plates




- Gamma sterilised with lid
- Microplate format for use with robots
- With notched corners : mark for positioning the plate
- Without notched corners : special robotics
- Applications :
  - Dot Blotting
  - Denaturation/hybridization/washing
  - Membrane storage

Cat. No.	Number of wells	Type of bottom	Specifications	Treatment	Units/ bag	Units/ carton	€ Excl. VAT/ carton
055048	1	Flat	Notched corners	TC	10	60	NC -
056199	1	Flat	-	TC	10	90	NC -

Nunc™ rectangular plates






- Sterile with cover
- Compatible with robotic systems
- ANSI standard format
- Ext. dimension L x W x H: 128 x 85.9 x 17.3 mm

Cat. No.	Number of wells	Type of bottom	Specifications	Treatment	Units/ bag	Units/ carton	€ Excl. VAT/ carton
056276	4	Flat	TC	10	100	60	NC -
056278	8	Flat	TC	10	100	90	NC -

Nunc™ transparent multi-well plates and microplates







- Sterilised by gamma rays
- In polystyrene
- Nunclon Delta Treated (TC)
- Nunclon VITA treated: for cultivation of so-called demanding cells such as human embryonic and pluripotent stem cells
- Minitray plates with cover: size L x W 84 x 59 mm, culture area per well 0.013 cm²
- Edge 2.0 plates: limits evaporation; evaporation rate ‹ 2 % on all 96 wells after 4 days of incubation
- 4 slots on the perimeter (each filled with sterile 0.5% agarose solution liquid), act as barriers to evaporation from the wells
- SBS format

Cat. No.	Number of wells	Type of bottom	Specifications	Treatment	Units/ bag	Units/ carton	€ Excl. VAT/ carton
055047	4	Flat	Lid	TC	4	120	NC -
055426	6	Flat	Lid	TC	1	75	NC -
055615	6	Flat	Lid	TC	5	85	NC -
056744	6	Flat	Lid	Nuclon Vita™	1	4	NC -
055428	12	Flat	Lid	TC	1	75	NC -
055429	24	Flat	Lid	TC	1	75	NC -
055430	24	Flat	Lid	TC	5	85	NC -
055431	48	Flat	Lid	TC	1	75	NC -
056204	48	Flat	Lid	TC	5	85	NC -
055435	60	Conical	Minitray	TC	10	150	NC -
055436	72	Conical	Minitray	TC	10	150	NC -
055260	96	Flat	Lid	TC	1	50	NC -
055262	96	Flat	Lid	TC	10	160	NC -
055432	96	Flat	Without lid	TC	1	50	NC -
055264	96	Flat	Without lid	TC	10	180	NC -
055434	96	Rounded	Lid	TC	1	50	NC -
055467	96	Rounded	Lid	TC	10	160	NC -
055433	96	Rounded	Without lid	TC	1	50	NC -
056891	96	Flat	Lid, Edge 2.0	TC	1	50	NC -
056892	96	Flat	Lid, Edge 2.0	TC	10	160	NC -
056893	96	Flat	Without lid, Edge 2.0	TC	1	50	NC -
056894	96	Flat	Without lid, Edge 2.0	TC	10	160	NC -

Sterile lid and adhesives

Cat. No.	Description	Pack unit.	Units/ carton	€ Excl. VAT/ carton
055610	Sterile polystyrene lid	1	60	NC -
055049	Sterile adhesive sheet	200	800	NC -
064031	Sterile adhesive sheet	1	50	NC -