CULTURE TUBES WITH DUAL POSITION CAP

Greiner Bio-One culture tubes with dual position cap








Volume (ml)	4.5		5	12	14		
Ø ext. (mm)	12.4	12.4	12	17	18	18	18
Ø int. (mm)	10	10	10	15	15	15	15
Height (mm)	75	75	75	76	95	95	95
Material	Polystyrene; HDPE cap		Polypropylene; LDPE cap	Polypropylene; LDPE cap	Polystyrene; HDPE cap		Polypropylene; LDPE cap
Sterile	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Bottom	Rounded	Rounded	Rounded	Rounded	Rounded	Rounded	Rounded
Caps	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Cap Type	Dual position snap						
Autoclavable	No	No	Yes (without cap)	No	No	No	Yes (without cap)
Max centrifugal speed (g)	5000	5000	34000	34000	5500	5500	34000
Graduation	No	No	Yes (1 ml increments)	Yes (1 ml to 10 ml increments)	No	No	Yes (1 ml increments)
Marking area	No	No	No	No	No	No	Yes
Colour	Translucent	Translucent	Translucent	Translucent	Transparent	Transparent	Translucent
Operating temperature	- 20 °C to +60 °C	- 20 °C to +60 °C	-196 °C to +121 °C	- 196 °C to +121 °C	- 20 °C to +60 °C	- 20 °C to +60 °C	-196 °C to +121 °C
Specifications	Free of RNase, DNase, human DNA and pyrogen; Non-cytotoxic	Culture-treated; Free of RNase, DNase, human DNA and pyrogen; Non-cytotoxic	Free of RNase, DNase, human DNA and pyrogen; Non-cytotoxic	Free of DNase et RNase, ADN humain and pyrogen; Non-cytotoxic	Culture-treated; Free of RNase, DNase, human DNA and pyrogen; Non-cytotoxic	Free of RNase, DNase, human DNA and pyrogen; Non-cytotoxic	Free of RNase, DNase, human DNA and pyrogen; Non-cytotoxic
Cat. No.	**120161**	**120160**	**115261**	**184261**	**-**	**-**	**187261**
Units/carton	1000 individually wrapped tubes						
€ Excl. VAT/carton	**NC -**	**NC -**	**NC -**	**NC -**	**-**	**-**	**NC -**
Cat. No.	**-**	**-**	**-**	**-**	**191160**	**191161**	**-**
Units/carton	750 individually wrapped tubes						
€ Excl. VAT/carton	**-**	**-**	**-**	**-**	**NC -**	**NC -**	**-**
Cat. No.	**120180**	**120190**	**115262**	**-**	**-**	**-**	**-**
Units/carton	80 bags of 25 tubes						
€ Excl. VAT/carton	**NC -**	**NC -**	**NC -**	**-**	**-**	**-**	**-**
Cat. No.	**-**	**-**	**-**	**-**	**-**	**191180**	**187262**
Units/carton	40 bags of 25 tubes						
€ Excl. VAT/carton	**-**	**-**	**-**	**-**	**-**	**NC -**	**NC -**

Nunc™ Culture tube with dual position cap

Volume (ml)	7
Ø (mm)	14
Height (mm)	100
Material	Polystyrene; polyethylene cap
Sterile	Yes
Bottom	Rounded
Cap Type	Snap; Dual position
Autoclavable	No
Max centrifugal speed (g)	-
Graduation	No
Marking area	No
Operating temperature	-
Specifications	Nunclon Delta® treatment
Cat. No.	**055053**
Units/carton	6 bags of 100 tubes
€ Excl. VAT/carton	**NC -**



Indicative prices. Visit www.dulis.nl