
MiniVertigel 2 vertical tank

- **Plate dimensions W x H : 10 x 10 cm**
- **Useful gel dimensions : 8 x 9 cm**
- **Total buffer volume : 250 ml**
- **Number of samples : 48**
- **Compatible with pre-cast gels in 10 x 10 cm cassettes**

Suggested generators
- PS 304 : 300 Volts - 400 mA - 120 Watts
- PS 305 : 300 Volts - 500 mA - 150 Watts with Timer and Gel Saver
- PS 503 : 500 Volts - 300 mA - 150 Watts with Timer and Gel Saver

The tank is supplied with :
- migration module ;
- lower tank ;
- cover with integral power cords ;
- 2 notched glass plates;
- 2 ceramic plates with notch ;
- 4 solid glass plates ;
- 1 blanking wall ;
- 2 combs 10 teeth 0.8 mm
- 4 spacers 0.8 mm ;
- 2 replacement seals ;
- casting base and seals.





Cat. No.	Description	€ Excl. VAT
049221	Complete MiniVertigel 2 for 2 gels 0,8 mm 10 x 10 cm with casting base	NC -
049241	Solid glass plate 10 x 10 cm	NC -
049242	Glass notched plate 10 x 10 cm	NC -
049243	Aluminium notched plate 10 x 10 cm	NC -


MidiVertigel 2 vertical tank

- **Plate dimensions W x H : 16 x 14 cm**
- **Useful gel dimensions : 14 x 12.5 cm**
- **Total buffer volume : 650 ml**
- **Number of samples : 48**

Suggested generators
- PS 305 : 300 Volts - 500 mA - 150 Watts with Timer and Gel Saver
- PS 503 : 500 Volts - 300 mA - 150 Watts with Timer and Gel Saver
- PS 608 : 600 Volts - 800 mA - 300 Watts with Timer and 3 displays

The tank is supplied with :
- migration module ;
- lower tank ;
- cover with integral power cords ;
- 4 notched glass plates;
- 4 solid glass plates ;
- 1 blanking wall ;
- 2 combs 10 teeth 1.5 mm
- 2 combs 15 teeth 1.5 mm
- 8 spacers 1.5 mm ;
- 2 replacement seals ;
- casting base and seals.





Cat. No.	Description	€ Excl. VAT
049222	MidiVertigel 2 for 2 gels 1.5 mm 16 x 14 cm with casting base	NC -
049263	Solid glass plate 16 x 14 cm	NC -
049264	Glass notched plate 16 x 14 cm	NC -


Mini Wide Vertigel 2 Vertical tank

- **Plate dimensions W x H : 20 x 10 cm**
- **Useful gel dimensions : 18 x 8 cm**
- **Total buffer volume : 650 ml**
- **Number of samples : 72**

Suggested generators
- PS 305 : 300 Volts - 500 mA - 150 Watts with Timer and Gel Saver
- PS 503 : 500 Volts - 300 mA - 150 Watts with Timer and Gel Saver
- PS 608 : 600 Volts - 800 mA - 300 Watts with Timer and 3 displays

The tank is supplied with :
- migration module ;
- lower tank ;
- cover with integral power cords ;
- 4 glass notched plates;
- 4 solid glass plates ;
- 1 blanking wall ;
- 2 combs 15 teeth 0.8 mm
- 2 combs 20 teeth 0.8 mm
- 8 spacers 0.8 mm ;
- 2 replacement seals ;
- casting base and seals.





Cat. No.	Description	€ Excl. VAT
049223	Mini Wide Vertigel 2 for 2 gels 0.8mm 20 x 10 cm with casting base	NC -
049286	Solid glass plate 20 x 10 cm	NC -
049287	Glass notched plate 20 x 10 cm	NC -



Find all compatible combs and accesories on www.dulis.nl

Indicative prices. Visit www.dulis.nl