
EasyPlate™Racks for Petri Dishes



- Increase the working area
- Eliminate the swinging effect of Petri dishes
- Good hold of the dishes thanks to an elastic system
- Supplied with marking labels
- Polycarbonate, Autoclavable

Two models available:
- for 48 (4 x 12) plastic Petri dishes 94 mm;
- for 84 (4 x 24) plastic Petri dishes 94 mm.


EasyPlate™ racks for 48 Petri dishes

Natural	Red	Yellow	Blue	€ Excl. VAT	Assorted	€ Excl. VAT/4 units
086645	086647	086648	086649	NC -	086651	NC -

EasyPlate™ racks for 84 Petri dishes

Natural	Blue	€ Excl. VAT
086646	086650	NC -


Rodac Racks for Petri dishes



- For 60 Petri dishes ø 65 mm or 42 Petri dishes ø 100 mm
- White polypropylene, with carrying handle
- Autoclavable at 121°C for 20 minutes
- Supplied flat, to be assembled
- L x W x H Cat.No. 065401: 26.6 x 16.5 x 17.8 cm
- L x W x H Cat.No. 088736: 33 x 21 x 17.8 cm

Cat. No.	Description	€ Excl. VAT
065401	Rack for 60 Petri dishes Ø 65 mm	NC -
088736	Rack for 42 Petri dishes Ø 100 mm	NC -

Incubation trays for Petri dishes



- Facilitate handling and incubation of Petri dishes
- White epoxy
- Dimensions : 349 x 134 x 22 mm
- Possible Petri dish diameters : 90, 60, 50 and 35 mm

Cat. No.	Description	€ Excl. VAT
085497	White epoxy incubation tray	NC -

Stackable incubation tray for Petri dishes




Polypropylene tray for 5 x 100 mm diameter Petri dishes or 15 x 60 mm diameter Petri dishes.

- Stackable trays to optimise incubator space
- L x W x H: 251 x 237 x 35 mm

Cat. No.	Description	€ Excl. VAT/ 3 units
086601	Stackable incubation tray for Petri dishes	NC -

Corning® Gosselin™Inoculating loops





Volume	1 µl	10 µl
Units/bag	20	20
Units/carton	9000	9000
Inoculating loop + ball	076835	076837
Inoculating loop + needle	076839	-
€ Excl. tax/carton	NC -	NC -

- Two models: 1 µl or 10 µl + Spherical extremity / or 1 µl or 10 µl + needle extremity (vaccinostyl)
- In soft polystyrene, softer with agar
- Ergonomic
- In re-sealable zip lock bag
- Sterile

Economic models

Cat. No.	Description	Pack unit.	Units/ carton	€ Excl. VAT
731165GO	Inoculating loop, 1 µl	20	3000	NC -
731175GO	Inoculating loop, 10 µl	20	3000	NC -