Cotton buds



- Wooden rod

Cat. No.	Number of heads	Rod size Ø x W (mm)	Ø head (mm)	Units/ box	€ Excl. VAT box
046027	2	2.2 x 80	6	100	NC -


Tongue depressor



- Wooden
- Smooth contact
- Rounded ends
- Non-sterile, single use

Cat.No	Description	Dimensions (mm)	Units/ box	€ Excl. VAT box
046009	Wooden tongue depressor	150 x 18	100	NC -
534801	Wooden tongue depressor ECO	150 x 18 x 2.2	100	NC -
534802	Wooden tongue depressor ECO single pack	180 x 21 x 2.35	100	NC -

Histobrush swabs




For endocervical sampling. Intact sampling Sterile brush swab for cell sampling, chlamydia testing, Papanicolaou collection, STD sampling, etc.

- Up to 70 times more effective than a traditional swab
- Sterile in individual packaging
- Length : 20 cm

Cat. No.	Description	€ Excl. VAT 100 units
252188	Histobruch sterile swabs	NC -

Nasopharyngeal swabs

Dry swabs



Rod	Material	ABS
	Length (mm)	153
	Ø (mm)	2.5
	Breakable (mm)	100
Tip	Material	Nylon
	Ø (mm)	3
Medium volume (ml)		3

- Rod material : ABS
- Total length : 150 mm
- Diameter : 2.5 mm
- Breaking height : 82.5 mm
- Flocked nylon tip
- Width : 3 mm
- Thickness : 3 mm
- Length : 16 mm
- Sterile
- CE and FDA
- Single pack

Cat. No.	Description	Units/ carton	€ Excl. VAT/ carton
011990	Swab in peelable bag	100	NC -

Swabs with medium

Cat. No.	Description	Units/ carton	€ Excl. VAT/ carton
010233	10 ml tube with 3 ml transport medium - preservation	500	NC -
010234	10 ml tube with 3 ml preservation medium - virus inactivation	500	NC -

Swab for DNA sampling



- Treated with ethylene oxide
- Tested for use in PCR
- Wooden rod, cotton tip
- Delivered with DNA analysis certificate

Cat.No	Description	Pack unit.	Units/ carton	€ Excl. VAT carton
020467A	Certified PCR DNA collection swabs	100	100	NC -

Swabs for fine sampling



- Laminated paper rod
- Mini cotton tip diameter 1.5 mm
- Gamma irradiation with Sal 10^{-6}

Cat. No.	Description	Units/ carton	€ Excl. VAT/ carton
039027B	Cotton swabs, paper rod, sterile tube	250	NC -