Acurex digital dispensers

- Fully protected mechanism inside the unit
- Graduated column retractable during refrigerated storage
- Autoclavable fully assembled
- Materials in contact with the liquid are inert
- Springless valve system
- Glass cylinder and piston with PTFE coating



Acurex 501 compact dispensers

4 models from 0.2 ml to 30 ml, supplied with yellow glass flask.

Cat. No.	Volume (ml)	Division (ml)	Capacity Flask (ml)	€ Excl. VAT
061027	0.2 - 2.0	0.1	250	NC -
061028	0.2 - 2.0	0.1	1000	NC -
061029	0.2 - 2.0	0.1	2000	NC -
061030	0.4 - 5.0	0.2	500	NC -
061031	0.4 - 5.0	0.2	1000	NC -
061032	0.4 - 5.0	0.2	2000	NC -
061033	1 - 10.0	0.2	1000	NC -
061034	1 - 10.0	0.2	2000	NC -
061035	1 - 30.0	1.0	2000	NC -

Accessories

Cat. No.	Description	€ Excl. VAT
061036	PTFE 60 cm extension tube for 1 and 2 litre tanks	NC -
061686	Replacement brown glass bottle, 2.5 litre	NC -

Eppendorf Dispensers Varispenser® 2 / Varispenser® 2X

- Dispensers for 0.5 to 100 ml flasks
- Good chemical resistance
- Suitable for a wide range of liquids, including corrosive and flammable liquids
- Varispencer 2X : return distribution valve prevents loss of reagent during priming as well as an easy calibration function
- Fully autoclavable
- Distribution block rotates 360 °
- Exit valve with secure ball closure that prevents liquid leakage when dispensing cannula is unassembled



				Varispenser® 2		Varispenser® 2X	
Volume (ml)	For neck (mm)	Reducers supplied	Increment (ml)	Cat.No	€ Excl. VAT	Cat.No	€ Excl. VAT
0.2 to 2	32	28, 40, 45	0.05	033580B	NC -	033586B	NC -
0.5 to 5	32	28, 40, 45	0.1	033581B	NC -	033587B	NC -
1 to 10	32	28, 40, 45	0.2	033592B	NC -	033588B	NC -
2.5 to 25	45	32, 38, 40	0.5	033610B	NC -	033589B	NC -
5 to 50	45	32, 38, 40	1	033584B	NC -	033590B	NC -
10 to 100	45	32, 38, 40	2	033585B	NC -	033591B	NC -

Accessory

Cat. No.	Description	€ Excl. VAT
999570	Threaded adapter for GL 45 to GL 32 bottles, ETFE	NC -

DispensMate reagent dispenser

Cat. No.	Volume (ml)	Graduation (ml)	Accuracy (%)	€ Excl. VAT
861070B	0.5 - 5	0.1	0.1 %	NC -
861071B	1 - 10	0.2	0.1 %	NC -
861073B	5 - 50	1.0	0.1 %	NC -

- Quick and precise volume control
- Drip cap for added safety
- Easy to clean
- Suitable for GL45 bottles (optional adapters)
- Fully autoclavable at 121°C



Accessory

Cat. No.	Description	€ Excl. VAT
861074	Bottle 1 l GL45 in amber glass	NC -

DispensMate Adapters

For bottle	S40	GL32	GL38	GL25	GL28
Cat. No.	861075	861076	861077	861078	861079B
€ Excl. VAT	NC -	NC -	NC -	NC -	NC -

Indicative prices. Visit www.dulis.nl