Finntips Tips


10 µl Flex


1200 µl Flex



- Very high precision tips
- Flex tips are more flexible, have better adhesion to the pipette and a sharper tip
- Natural colour, colour-coded racks
- Sterile Finntip certified free of endotoxins, DNA, DNase and RNase (standard, filter and Biocon tips)
- Ext. model : extra long
- Racks with cover

Designation	Volume	Lenght (cm)	Code colour		Bag	on rack	on rack stérile	Refill		Box	Bulk eco	BioCon single sterile
Finntip 10 µl	0,2 - 10 µl	3,2	Pink	Pack.	1,000 units	10 x 96	10 x 96	10 x 192*	20 x 192***	-	-	100
				Cat. No.	047042	047043	047044	047260	047257	-	-	047157
				€ Excl. VAT	NC -	NC -	NC -	NC -	NC -	-	-	NC -
Finntip 10 µl Flex	0,2 - 10 µl	3,3	Pink	Pack.	1,000 units	10 x 96	10 x 96	20 x 96**	20 x 96	-	-	-
				Cat. No.	047950A	047946	047947	047948	047949	-	-	-
				€ Excl. VAT	NC -	NC -	NC -	NC -	NC -	-	-	-
Finntip 20 µl	0.2 - 20 µl	3.3	Purple	Pack.	1,000 units	10 x 384	10 x 384	-		-	-	-
				Cat. No.	047045	047046	047047	-		-	-	-
				€ Excl. VAT	NC -	NC -	NC -	-		-	-	-
Finntip 50 µl (For Finnpipette 10 µl and 16 channel 50 µl)	0,2 - 50 µl	4,8	Turquoise	Pack.	1,000 units	10 x 384	10 x 384	-		-	-	-
				Cat. No.	047048	047602	047726	-		-	-	-
				€ Excl. VAT	NC -	NC -	NC -	-		-	-	-
Finntip 200 µl ext.	5 - 200 µl	7.8	Orange	Pack.	-	10 x 96	10 x 96	-		Box of 400	-	-
				Cat. No.	-	047052	047053	-		047051	-	-
				€ Excl. VAT	-	NC -	NC -	-		NC -	-	-
Finntip 200 µl Flex	5 - 200 µl	6	Yellow	Pack.	1,000 units	10 x 96	10 x 96	20 x 96**	20 x 96	-	-	-
				Cat. No.	047274	047275	047276	047277	047278	-	-	-
				€ Excl. VAT	NC -	NC -	NC -	NC -	NC -	-	-	-
Finntip 250 µl Universal	0.5 - 250 µl	5.2	Yellow	Pack.	1,000 units	10 x 96	10 x 96	20 x 96	10 x 96*	box of 500	20 000	100
				Cat. No.	047055	047057	047058	047059	047245	047054	047056	047158
				€ Excl. VAT	NC -	NC -	NC -	NC -	NC -	NC -	NC -	NC -
Finntip 250 µl wide (large opening)	0.5 - 250 µl	5.2	White	Pack.	1,000 units	10 x 96	10 x 96	-		-	-	-
				Cat. No.	047758	047760	047761	-		-	-	-
				€ Excl. VAT	NC -	NC -	NC -	-		-	-	-
Finntip 300 µl	5 - 300 µl	5.2	Orange	Pack.	1,000 units	10 x 96	10 x 96	10 x 192		-	20 000	-
				Cat. No.	047060	047062	047063	047064		-	047061	-
				€ Excl. VAT	NC -	NC -	NC -	NC -		-	NC -	-
Finntip Flex 300 µl	5 - 300 µl	6.0	Orange	Pack.	1,000 units	10 x 96	10 x 96	20 x 96	20 x 96**	-	-	-
				Cat. No.	047755	047727	047728	047730	047729	-	-	-
				€ Excl. VAT	NC -	NC -	NC -	NC -	NC -	-	-	-
Finntip 1,000 µl	100 - 1,000 µl	7.1	Blue	Pack.	1,000 units	10 x 96	10 x 96	5 x 192		box of 200	10 000	50
				Cat. No.	047065	047247	047248	047249		047067	047066	047159
				€ Excl. VAT	NC -	NC -	NC -	NC -		NC -	NC -	NC -
Finntip 1,000 µl ext.	100 - 1,000 µl	10.5	Blue	Pack.	1,000	The 5 x 96	The 5 x 96	-		-	-	-
				Cat. No.	047071	047250	047251	-		-	-	-
				€ Excl. VAT	NC -	NC -	NC -	-		-	-	-
Finntip 1,000 µl wide (large opening)	100 - 1,000 µl	7,2	White	Pack.	400 units	10 x 96	10 x 96	-		-	6500	-
				Cat. No.	047762	047764	047765	-		-	047763	-
				€ Excl. VAT	NC -	NC -	NC -	-		-	NC -	-
Finntip 1,000 µl Flex	100 - 1,000 µl	8.7	Blue	Pack.	1,000	10 x 96	10 x 96	2 x 8 x 96	2 x 8 x 96****	-	-	-
				Cat. No.	047849	047848	047852	047850	047851	-	-	-
				€ Excl. VAT	NC -	NC -	NC -	NC -	NC -	-	-	-
Finntip 1200 µl Flex	50 - 1200 µl	8,7	Turquoise	Pack.	1,000	10 x 96	10 x 96	2 x 8 x 96	2 x 8 x 96****	-	-	-
				Cat. No.	047858	047756	047757	047953	047952	-		-
				€ Excl. VAT	NC -	NC -	NC -	NC -	NC -	-	-	-
Finntip 5 ml	1 - 5 ml	14,7	Green	Pack.	500 units	5 x 54	5 x 54	-		Box with 75	3000	-
				Cat. No.	047075	047100	047101	-		047074	047076	-
				€ Excl. VAT	NC -	NC -	NC -	-		NC -	NC -	-
Finntip 10 ml	2 - 10 ml	15	Red	Pack.	100 units	5 x 24	5 x 24	-		40 units	-	-
				Cat. No.	047080	047081	047082	-		047079	-	-
				€ Excl. VAT	NC -	NC -	NC -	-		NC -	-	-
Finntip 10 ml Flex Long	1 - 10 ml	26,7	Red	Pack.	100 units	-	50 (bulk)	-		-	-	-
				Cat. No.	047911	-	050716	-		-	-	-
				€ Excl. VAT	NC -	-	NC -	-		-	-	-
Finntip Multistepper	1500 µl	8,5	Turquoise	Pack.	400 units	10 x 96	10 x 96	-		-	-	-
				Cat. No.	047170	047171	047172	-		-	-	-
				€ Excl. VAT	NC -	NC -	NC -	-		-	-	-
Finntip PDP 25	0,2 - 25 µl	3,7	Blue light	Pack.	960 units	10 x 96	-	-		-	-	-
				Cat. No.	047139	047140	-	-		-	-	-
				€ Excl. VAT	NC -	NC -	-	-		-	-	-

Tips not shown to scale. Refer to the lengths (in cm) in the table.
* With 4 racks (Starter Kit) ** 20 refills of 96 + 1 rack of 96 *** Loaded tray (without rack) **** 8 trays of 96 x 2 + 1 full rack