ClearLine® Xtip4 Tips for Rainin LTS pipettes

- Graduated polypropylene tips
- Compatible with the LTS Micropipettes
- Flexibility of the tip for lower insertion force and good sealing
- Locking device on the tip to reduce the ejection force
- No droplet formation at the tip
- Low binding



Pyro-genic | DNase RNase | Nucleic acid | Endo-toxins | PCR-inhibitor

Xtip4 tips without filter

		On a rack				RechargeTipEject	
		Non-sterile		Sterile		Sterile	
Vol. (µl)	Units/ carton	Cat. No.	€ Excl. VAT/ carton	Cat. No.	€ Excl. VAT/ carton	Cat. No.	€ Excl. VAT/ carton
20	10 x 96	713184B	NC -	713181B	NC -	713193B	NC -
250	10 x 96	713185B	NC -	713182B	NC -	713194B	NC -
300	10 x 96	713230	NC -	713228	NC -	-	-
1000	8 x 96	713186B	NC -	713183B	NC -	713195B	NC -
1200	8 x 96	713231	NC -	713229	NC -	713233	NC -

Xtip4 sterile filter tips

Volume (µl)	Units/ carton	Cat. No.	€ Excl. VAT/ carton
20	10 x 96	713178B	NC -
200	10 x 96	713179B	NC -
300	10 x 96	713221	NC -
1000	8 x 96	713180B	NC -
1200	8 x 96	713222	NC -

Accessories

Cat. No.	Description	Units/ carton	€ Excl. VAT/ carton
713196B	Empty box for refill Xtip4 20 - 300 µl	50	NC -
713197B	Empty box for refill Xtip4 1000 - 1200 µl	40	NC -

2 in 1 system!


Work with theClearLine® non filter and filter range and reduce your waste : consider reusing your hinged racks with the refill systems

GoGreen

Instructions for use : Tip Eject refill

- Tower containing 10 trays of 96 sterile tips
- Easy one-handed operation
- Recyclable packaging

- 1°/remove the protective rack at the base of the column
- 2°/place the tower on the empty rack
- 3°/press the top - a “click” indicates that the tray is correctly positioned in the rack
- 4°/replace the column on the protective rack



Indicative prices. Visit www.dulis.nl