Sterile Tips, unit packaging

Gamma sterilised micropipette tips, packed in individual bags.

- **Ideal for small pipetting runs**
- **Risks of contamination completely eliminated**
- **No autoclaving**
- **Do not take up space under the hood**
- **Safety and ease of use**



Standard Tips

Cat. No.	Tip	Volume (µl)	Compatibility	Units/ box	€ Excl. VAT/ box
077209	Standard yellow	200	Universal	500	**NC -**
077169	Standard blue	1000	Universal	500	**NC -**

High precision tips

Free of heavy metals.

DNA · Pyro-genic · Endo-toxins · DNase RNase

Cat. No.	Tip	Volume (µl)	Compatibility	Units/ box	€ Excl. VAT/ box
949472	Yellow	200	Universal	200	**NC -**
949473	Natural graduated	200	Universal	200	**NC -**
949474	Natural extra-long	1250	Universal	200	**NC -**

5 ml and 10 ml tips

5 ml and 10 ml tips for Gilson, Eppendorf and Labsystems.

- **In bulk or in refillable rack with lid**
- **Excellent quality**
- **With or without filter**



		Without filter - bulk			Without filter - on rack			Without filter - on sterile rack			With Filter - on sterile rack		
Volume	Length (mm)	Cat.No	Units/ bag	€ Excl. VAT/ bag	Cat.No	Units/ bag	€ Excl. VAT/ bag	Cat.No	Units/ bag	€ Excl. VAT/ bag	Cat.No	Units/ bag	€ Excl. VAT/ bag
5 ml for Gilson	125	**045413**	250	**NC -**	**045414**	10 x 50	**NC -**	**045415**	10 x 50	**NC -**	**045561**	10 x 50	**NC -**
5 ml for Eppendorf, Sartorius	155	**045530**	250	**NC -**	**045531**	10 x 50	**NC -**	**045532**	10 x 50	**NC -**	**045562**	10 x 50	**NC -**
10 ml for Socorex	165	**045564**	200	**NC -**	**045565**	10 x 25	**NC -**	**048060**	10 x 25	**NC -**	**-**	-	**-**
10 ml for Gilson, Rainin Eppendorf	161	**048026**	200	**NC -**	**048027**	10 x 25	**NC -**	**048028**	10 x 25	**NC -**	**048029**	10 x 25	**NC -**

Tips for Pipet-Lite Rainin

	20 µl		200 µl		1000 µl	
	Non-sterile	Sterile	Non-sterile	Sterile	Non-sterile	Sterile
Standard without filter						
Cat. No.	**-**	**038373**	**038375**	**038377**	**038380**	**038381**
Units/box	-	960	960	960	800	800
€ Excl. VAT/box	**-**	**NC -**	**NC -**	**NC -**	**NC -**	**NC -**
Standard with ART filter						
Cat. No.	**-**	**038374**	**-**	**038379**	**-**	**038382**
Units/box	-	960	-	960	-	800
€ Excl. VAT/box	**-**	**NC -**	**-**	**NC -**	**-**	**NC -**
Low retention without filter						
Cat. No.	**038383**	**038384**	**038386**	**038387**	**-**	**-**
Units/box	960	960	960	960	-	-
€ Excl. VAT/box	**NC -**	**NC -**	**NC -**	**NC -**	**-**	**-**
Low retention with ART filter						
Cat. No.	**-**	**038385**	**-**	**038388**	**-**	**038391**
Units/box	-	960	-	960	-	800
€ Excl. VAT/box	**-**	**NC -**	**-**	**NC -**	**-**	**NC -**



- **Specially designed tips for the Rainin Pipet-Lite pipette**
- **Very low adsorption thin walls**
- **On rack, with or without ART hydrophobic filter**
- **Truncated cone tip for effortless attachment and ejection of the tip**
- **Standard or Low-retention with low surface adhesion**

Indicative prices. Visit www.dulis.nl