ESTY (Broth)

For the isolation and enumeration of lactic streptococci in yoghurt and dairy products.

➤ **Storage : 2 - 25 °C**

➤ **Incubation : 35 ± 2 °C for 24 - 48 h**

Cat. No.	Description	Pack.	€ Excl. VAT
777247	Dehydrated medium for ESTY broth - Condalab	500 g	NC -

ESTY (agar)

For the enumeration of thermophilus streptococci in yoghurt.

➤ **Storage : 2 - 25 °C**

➤ **Incubation : 35 ± 2 °C for 24 - 48 h**

Cat. No.	Description	Pack.	€ Excl. VAT
777248	Dehydrated medium for ESTY agar - Condalab	500 g	NC -

Eugon LT100 (broth)

For the enrichment of micro-organisms. Recommended for the enumeration of aerobic mesophilic bacteria, yeast and moulds and Pseudomonas aeruginosa, Staphylococcus aureus.

➤ **Storage : 2 - 25 °C (Dehydrated medium and ready-to-use medium)**

➤ **Incubation : 30 - 35 °C for 24 to 72 h**

Dehydrated medium

Cat. No.	Description	Pack.	€ Excl. VAT
778403	Dehydrated medium for eugon broth LT100 - Condalab	500 g	NC -

Ready-to-use media

Cat. No.	Description	Pack.	€ Excl. VAT
693300	9 ml tube of eugon broth LT100 - ClearLine®	100	NC -
693401	Eugon Broth 90 ml bottle LT100 - ClearLine®	10	NC -

Ewing malonate (broth)

For the differentiation of coliforms.

➤ **Storage : 2 - 25 °C**

➤ **Incubation : 35 ± 2 °C for 18 - 48 h**

Cat. No.	Description	Pack.	€ Excl. VAT
777251	Dehydrated medium for modified Ewing malonate broth - Condalab	500 g	NC -

Yeast extract

An ingredient used as a source of nitrogen in media.

➤ **Storage : 2 - 25 °C**

Cat. No.	Description	Pack.	€ Excl. VAT
777474	Yeast extract - Condalab	500 g	NC -

Malt extract

An ingredient used as a source of nitrogen in media.

➤ **Storage : 10 - 35 °C**

Cat. No.	Description	Pack.	€ Excl. VAT
777462	Malt extract - Condalab	500 g	NC -

Fraser

For the detection and enumeration of Listeria in food products.

➤ **Storage : 2 - 25 °C (for the base) and 2 - 8 °C (for the ready-to-use medium)**

➤ **Incubation : 30°C ± 1°C for 24 hours ± 2 hours**

Dehydrated medium

Cat. No.	Description	Pack.	€ Excl. VAT
777547	Fraser Broth Dehydrated Medium - ISO 11133 / ISO 11290 compliant - Condalab	500 g	NC -

Supplement

➤ **Ferric ammonium citrate. To be reconstituted in 5 ml of warm sterile distilled water**

Cat. No.	Description	Pack.	€ Excl. VAT
777782	Ferric ammonium citrate supplement - Condalab	10 bottles	NC -

Ready-to-use media

Cat. No.	Description	Pack.	€ Excl. VAT
693409	10 ml tube of Fraser broth in accordance with ISO 11290-1 - ClearLine®	100	NC -

Half Fraser

For the detection and enumeration of Listeria in food products.

➤ **Storage : 2 - 25 °C (for the base) and 2 - 8 °C (for the ready-to-use medium)**

➤ **Incubation (according to standard 11290) : 30 °C for 25 hours ± 1 hour**

Dehydrated medium

Cat. No.	Description	Pack.	€ Excl. VAT
777852	Fraser demi broth dehydrated medium - ISO 11133 / ISO 11290 compliant - Condalab	500 g	NC -

Supplement

➤ **Ferric ammonium citrate. To be reconstituted in 5 ml of warm sterile distilled water**

Cat. No.	Description	Pack.	€ Excl. VAT
777782	Ferric ammonium citrate supplement - Condalab	10 bottles	NC -

Ready-to-use media

Cat. No.	Description	Pack.	€ Excl. VAT
693460	225 ml bottle of Half Fraser broth - ClearLine®	6	NC -
993950B	3 litres pouch of Half Fraser broth - ClearLine®	4	NC -
993956B	5 litres pouch of Half Fraser broth - ClearLine®	2	NC -
994512	5 litres pouch of Half Fraser broth - ClearLine®	55	NC -

Lactose gelatin



For confirmation of clostridium perfringens.

➤ **Storage : 2 - 25 °C**

➤ **Incubation : 35 ± 2 °C for 18 - 48 h**

Cat. No.	Description	Pack.	€ Excl. VAT
777256	Dehydrated medium gelatine lactose broth - Condalab	500 g	NC -