Micrewtube® without cap Simport®

	Non-graduated									
	Standard					For insert				
										
Volume (ml)	0.5	1.5		2		0.5	1.5		2	
Diameter (mm)	11									
Height (mm)	44									
Material	Polypropylene									
Sterile	No									
Autoclavable	Yes									
Max centrifugal speed (g)	17000	17000	20000	17000	20000	17000	17000	20000	17000	20000
Skirt	Yes	Yes	No	Yes	No	Yes	Yes	No	Yes	No
Graduation	No									
Marking area	No									
Operating temperature	-196 °C to +121 °C									
Specifications	Free of release agents, plasticizers and biocides; free of DNases, DNA RNases and pyrogens					Very low adhesion material for maximum sample recovery; free of release agents, plasticizers and biocides; free of DNases, DNA RNases and pyrogens				
Units/box	1000 microtubes									
Ref. microtube only	**039171**	**039173**	**039172**	**039175**	**039174**	**039181**	**039183**	**039182**	**039185**	**039184**
€ Excl. VAT/box	**NC -**	**NC -**	**NC -**	**NC -**	**NC -**	**NC -**	**NC -**	**NC -**	**NC -**	**NC -**

	Graduated, low adhesion					Amber				
										
Volume (ml)	0.5	1.5		2		0.5	1.5		2	
Diameter (mm)	11									
Height (mm)	44									
Material	Polypropylene									
Sterile	No									
Autoclavable	Yes									
Max centrifugal speed (g)	17000	17000	20000	17000	20000	17000	17000	20000	17000	20000
Skirt	Yes	Yes	No	Yes	No	Yes	Yes	No	Yes	No
Graduation	Yes	Yes	Yes	Yes	Yes	No	No	No	No	No
Marking area	Yes	Yes	Yes	Yes	Yes	No	No	No	No	No
Operating temperature	-196 °C to +121 °C									
Specifications	Free of release agents, plasticizers and biocides; free of DNases, DNA RNases and pyrogens.					Free of release agents, plasticizers and biocides; Free of DNases, DNA RNases and pyrogens; For light-sensitive samples.				
Units/box	1000 microtubes									
Ref. microtube only	**039176**	**039178**	**039177**	**039180**	**039179**	**-**	**-**	**-**	**-**	**-**
Cat. No. single amber microtube	**-**	**-**	**-**	**-**	**-**	**039186**	**039188**	**039187**	**039190**	**039189**
€ Excl. VAT/box	**NC -**	**NC -**	**NC -**	**NC -**	**NC -**	**NC -**	**NC -**	**NC -**	**NC -**	**NC -**

Indicative prices. Visit www.dulis.nl