CAPS

TIP

To choose the right cap: measure the internal diameter of the tube for inverted caps and the external diameter for capping caps

Caps for 12 and 17 mm tubes

Simport



For tubes Ø (mm)	12	17
Units/box	Box of 4000 caps	4 bags of 1000 caps
White	**039974**	**039398**
Yellow	**039975**	**039397**
Red	**039973**	**039395**
Green	**039972**	**039394**
Blue	**039971**	**039396**
Natural	**039970**	-
€ Excl. VAT/box	**NC -**	**NC -**

Re-entrant caps

Thermo SCIENTIFIC | nunc



For tubes Ø (mm)	11	12	15
Units/box	6 bags with 600 Plugs	6 bags with 500 Plugs	8 bags with 300 Plugs
Natural	**056045**	**056047**	**056050**
€ Excl. VAT/box	**NC -**	**NC -**	**NC -**

Double position cap

labcon



For tubes Ø (mm)	12	17
Units/box	3 bags of 1000 caps	Bag of 1000 caps
Natural	**949919**	**949958**
€ Excl. VAT/box	**NC -**	**NC -**

Re-entrant caps

labcon



For tubes Ø (mm)	17
Units/box	2 bags with 500 Plugs
Natural	**949877**
€ Excl. VAT/box	**NC -**

Caps



For tubes Ø (mm)	outer Ø 16 mm and inner Ø 13.5 mm
Packaging unit	500
Units/box	12
White	**483041**
€ Excl. VAT/box	**NC -**

Plug cap

labcon



For tubes Ø (mm)	12	
Units/box	3 bags of 1000 caps	6 bags of 500 caps
White	**949921**	-
Yellow	**949922**	-
Red	**949923**	-
Green	**949926**	-
Blue	**949925**	-
Natural	**949920**	-
Assorted	-	**949927**
€ Excl. VAT/box	**NC -**	**NC -**