Expanded polystyrene racks



- For the transport of 50 x 12 mm diameter tubes
- Dimensions : 146 x 73 x 40 mm
- A lightweight cardboard rack with dimensions of 155 x 155 x 45 mm enables two polystyrene boxes to be gathered together


Expanded polystyrene racks

White	Blue	Green	Red	€ Excl. VAT
049725	049726	049727	049728	NC -

Cardboard rack

Cat. No.	Description	€ Excl. VAT
049729	Cardboard Rack	NC -

Stainless steel cryobox with lid



Cat. No.	Dimensions with lid W x D x H (mm)	€ Excl. VAT
140076	133 x 133 x 50	NC -
140077	133 x 133 x 75	NC -
140078	133 x 133 x 100	NC -

Cryotube clamp



- For cryotubes
- Made of polypropylene

Cat. No.	Description	€ Excl. VAT/ 12 units
029111	Clamp for cryotubes	NC -

Supporting rod for cryo-preservers and Cryoflex sheath

Cryoflex sheath

Sheath for holding cryotubes, to be welded to prevent accidental entry of liquid nitrogen into the tube during total immersion. Length : 50 cm.



Cat. No.	Description	€ Excl. VAT 15 sheaths
055029	Cryoflex sheath	NC -

Thermo SCIENTIFIC nunc



Aluminium rod to hold 3 to 6 cryotubes depending on their height. Length : 28 cm.

Cat. No.	Description	€ Excl. VAT 50 units
055030	Rod holder	NC -

Labels and label holder



Transparent polypropylene label holder with one adhesive side.
To be placed on one side of the freezer box when it is completely dry and at room temperature.

- For temperatures down to -90 °C by stapling the label holder
- Adhesive side works between -10 °C and +70 °C

Cat. No.	For	Units/ box	€ Excl. VAT box
049616	Label holder	20	NC -
049617	Labels	5 x 20	NC -