Polypropylene bottles, wide neck

Azlon® wide neck polypropylene bottles



- Very good sealing
- Polypropylene bottles, autoclavable at + 121 °C
- Good chemical resistance

DWK LIFE SCIENCES

		10 units		BULK version		
Vol. (ml)	Ø x H (mm)	Cat. No.	€	Cat. No.	Units/ box	€/ box
30	32 x 70	BWP0030N	NC -	785409	1200	NC -
60	40 x 82	BWP0060N	NC -	761866	740	NC -
100	51 x 89	232661	NC -	-	-	-
125	51 x 104	BLP0125P	NC -	761867	350	NC -
150	51 x 116	232662	NC -	761868	350	NC -
250	60 x 140	BWP0250P	NC -	761869	200	NC -
500	75 x 161	BWP0500P	NC -	761872	100	NC -
1000	90 x 205	BWP1000P**	NC -	761870	60	NC -

** 5 units

Azlon® Amber polypropylene bottles with wide neck



- Ideal for photosensitive products
- Opaque
- Very good sealing
- Polypropylene bottles, autoclavable at + 121 °C

DWK LIFE SCIENCES

Cat. No.	Volume (ml)	ø (mm)	Height (mm)	Units/ box	€/ box
232804	60	40	82	10	NC -
723866	150	51	116	10	NC -
232805	250	60	140	10	NC -
232545	500	75	161	10	NC -
232806	1000	90	205	5	NC -
232761	2000	120	242	5	NC -

Polypropylene standardised wide neck graduated flasks



WEB

Enter the Cat. No. to find our autoclaving procedure

- With graduations
- Durable and long-lasting, autoclavable
- Certified for food contact
- Neck and thread according to DIN 13316 and 168

Vol. (ml)	Grad. (ml)	ø x H (mm)	Neck ø (mm)	DIN	Units/ box	Cat. No.	€/ box
50	10	38 x 88	24	GL32	100	391273	NC -
100	20	48 x 105	24	GL32	100	391274	NC -
250	25	60 x 140	38	GL45	50	391275	NC -
500	100	75 x 170	38	GL45	25	391276	NC -
1000	100	95 x 206	55	GL63	17	391277	NC -
2000	100	120 x 252	55	GL63	15	391278	NC -

WEB

Find VITgrip® screw bottles on dulis.nl

Nalgene™large round bottles with wide opening



In Autoclavable polypropylene. For powder transport or water sampling.

Thermo SCIENTIFIC NALGENE

Capacity (l)	Ø x H (mm)	cap Ø (mm)	Cat. No.	€
2	119 x 244	100	029032A	NC -
4	153 x 288	100	029534A	NC -

Nalgene™High-strength bottles



- Thick walls, reinforced closure (watertightness guaranteed)
- For the safe storage of media, biological or chemical liquids, dry chemicals
- Excellent chemical resistance, ideal for vacuum applications

Thermo SCIENTIFIC NALGENE

Volume (ml)	250	1000	2000	4000	5000
ø x H (mm)	74 x 135	93 x 226	119 x 260	155 x 338	156 x 414
Capuchon type	53B	53B	53B	83B	83B
Cat. No.	029307	029305	029306	029296	028162
Units/bag	1	6	2	1	1
€/bag	NC -	NC -	NC -	NC -	NC -