Sampling pots with screw cap, labelled, made of polystyrene

Sampling pots with metal lid

- In polystyrene, with metal cap and high sterility seal
- CE marking in accordance with the Directive on Medical Devices for in vitro diagnosis (98/79/EC)
- Aseptic
- 60 ml pots guaranteed 95 kPa



		With printed inscription surface			With blank inscription surface		
Volume (ml)	Ø x H (mm)	Cat. No.	Units/ box	€/ box	Cat. No.	Units/ box	€/ box
60	44 x 61	**761606**	300	**NC -**	**761607**	300	**NC -**
100	50 x 76	**761609**	200	**NC -**	**761610**	200	**NC -**
150	55 x 107.5	**761611**	120	**NC -**	**761612**	120	**NC -**
250	65 x 121	**761613**	50	**NC -**	**761614**	50	**NC -**
250	65 x 121	**761617***	50	**NC -**	-	-	-

* On a tray

Aseptic Bijou bottle 7 ml with screw cap

Flat-bottomed polystyrene tube with screw cap.

- Volume : 7 ml
- Aseptic
- ø x H : 17 x 50 mm
- White screw cap



Cat. No.	Description	Units/ box	€/ box
080171	Tube with label	700	**NC -**

Wheaton® Universal flasks

- Watertight closure.
- Dimensions ø x H : 30.5 x 93.5 mm
- Volume : 30 ml
- Sterile

DWK LIFE SCIENCES



Cat. No.	Material	Label	Boric acid	Packaging	€
055730	Polystyrene	Virgin	No	8 bags of 50	**NC -**
055732	Polystyrene	Printed	No	8 bags of 50	**NC -**

Universal skirted containers 30 ml Sterilin

Label	Label marking	Sterile	Units/ box	Cat. No.	€/ box
With	Yes	No	400	**761618**	**NC -**
With	No	No	400	**761619**	**NC -**
With	Yes	Yes	400	**080005E**	**NC -**

* Cat. No. sold with a blue plastic spoon

- With polypropylene screw caps and high sterility seal
- Conical skirted bottom
- Can be centrifuged at 3889 g
- 30 ml tubes guaranteed 95 kPa
- ø x H : 24.5 x 89 mm
- CE marking in accordance with the Directive on Medical Devices for In Vitro Diagnostics (98/79/EC)



All prices are indicative and exclude VAT. Visit www.dulis.nl