## Wheaton® Straight Edge Jars

- With white polypropylene cap, PTFE-faced foamed polyethylene seal
- For sampling dry goods, powders
- Made of soda lime glass, USP type III compliant



| Cat. No. | Capacity (ml) | ø x H (mm) | Cap | Units/ box | €/ box |
|---|---|---|---|---|---|
| 048715 | 60 | 55 x 48 | 53-400 | 24 | NC - |
| 048716 | 125 | 60 x 68 | 58-400 | 24 | NC - |
| 048717 | 250 | 73 x 88 | 70-400 | 12 | NC - |
| 048718 | 500 | 91 x 95 | 89-400 | 12 | NC - |
| 048719 | 1000 | 95 x 170 | 89-400 | 12 | NC - |

## Wheaton®Square Bottles

Space-saving bottles. Made of soda lime glass, USP type III compliant. With white polypropylene cap and PTFE-faced foamed polyethylene seal.



| Cat. No. | Volume (ml) | ø x H (mm) | Cap | Units/ box | €/ box |
|---|---|---|---|---|---|
| 048720 | 15 | 26 x 62 | 20-400 | 48 | NC - |
| 048721 | 30 | 31 x 72 | 24-400 | 48 | NC - |
| 048722 | 60 | 39 x 87 | 28-400 | 48 | NC - |
| 048723 | 125 | 45 x 111 | 33-400 | 24 | NC - |
| 048724 | 250 | 56 x 137 | 43-400 | 24 | NC - |
| 048725 | 500 | 68 x 167 | 48-400 | 24 | NC - |

## Wheaton® Graduated Storage Bottles with cap




- Made of soda lime glass, USP type III compliant
- For dry products and powders, double graduation in Oz and ml

| Cat. No. | Volume (ml) | ø x H (mm) | Cap | Units/ box | €/ box |
|---|---|---|---|---|---|
| With white polypropylene cap and PTFE-faced foamed polyethylene seal. | | | | | |
| 048695 | 30 | 34 x 68 | 33-400 | 48 | NC - |
| 048696 | 60 | 42 x 83 | 38-400 | 48 | NC - |
| 048697 | 125 | 51 x 102 | 48-400 | 24 | NC - |
| 048698 | 250 | 62 x 127 | 58-400 | 24 | NC - |
| 048699 | 500 | 76 x 145 | 70-400 | 24 | NC - |
| With black phenolic cap and autoclavable rubber seal | | | | | |
| 048700 | 30 | 34 x 68 | 33-400 | 48 | NC - |
| 048701 | 60 | 42 x 83 | 38-400 | 48 | NC - |
| 048702 | 125 | 51 x 102 | 48-400 | 24 | NC - |
| 048703 | 250 | 62 x 127 | 58-400 | 24 | NC - |
| Without cap | | | | | |
| 048704 | 30 | 34 x 68 | 33-400 | 432 | NC - |
| 048705 | 60 | 42 x 83 | 38-400 | 288 | NC - |
| 048706 | 125 | 51 x 102 | 48-400 | 144 | NC - |
| 048707 | 250 | 62 x 127 | 58-400 | 96 | NC - |
| 048708 | 500 | 76 x 145 | 70-400 | 48 | NC - |

## Amber Glass Powder Bottles

- Amber Soda-lime glass
- Supplied without cap

| Cat. No. | Volume (ml) | ø x H (mm) without cap | Cap | Units/ box | €/ box |
|---|---|---|---|---|---|
| 674481 | 30 | 38.10 x 57.90 | 33-400 | 127 | NC - |
| 674413 | 60 | 44 x 75.40 | 38-400 | 72 | NC - |
| 674414 | 100 | 48.40 x 88.10 | 38-400 | 72 | NC - |
| 674415 | 125 | 52.80 x 94.50 | 38-400 | 68 | NC - |
| 674419B | 500 | 81.10 x 148.70 | 58-400 | 28 | NC - |

### Polypropylene caps only

| Cat. No. | Description | € 100 units |
|---|---|---|
| 674482 | Black cap 33-400 | NC - |
| 674483 | Black cap 38-400 | NC - |
| 674416 | White cap 38-400 | NC - |
| 674418 | Black cap 48-400 | NC - |
| 674420 | Black cap 58-400 | NC - |


(Supplied without cap)



All prices are indicative and exclude VAT. Visit www.dulis.nl