Anios RN Neutralizer



Machine rinser and neutralizer. Neutralization of alkaline residues from detergents and water during the automatic machine washing cycle of medical-surgical devices and instruments.

- **Suitable for short tank washing cycles**
- **Complexing effect of minerals**
- **Traceable by conductimetry (1.77 mS/cm)**
- **Formulated with phosphoric acid (pH 2.5)**
- **Directions for use: dilute between 0.1 and 0.5% (1 to 5 ml/l); rinse with water at the end of the cycle**

Cat. No.	Description	€
059941A	ANIOS RN - 5 l canister	NC -
059942	ANIOS RN - 10 l canister	NC -

Anios DLH Detergent



Instrumentation detergent. Cleaning of devices and equipment.

- **Alkaline products**
- **Use with Anios RN**
- **Can be used in any water quality**
- **Non-corrosive**
- **pH 13.1**
- **Stability up to 95°C**
- **Instructions for use: dilution between 0.3 and 0.5%; optimum temperature between 55°C and 70°C**

Cat. No.	Description	€
059948A	ANIOS DLH - 5 l canister	NC -

Specific detergent for glassware Lancer LLL



- **Machine washable up to 95°C**
- **Applicable to all types of glassware: laboratory, foodstuffs, baby bottles**
- **Sodium hydroxide based, environmentally friendly in normal use**
- **In 10 litre cans**

NEW FORMULA :

- **Phosphate-free**
- **Non-toxic**

Cat. No.	Description	€
331700B	Liquid detergent LLL 10 litres	NC -

Lancer NLL liquid neutralizer



- **Specific rinsing solution for glassware, allows to obtain a neutral pH after washing**
- **Applicable to all types of glassware: laboratory, foodstuffs, baby bottles**
- **Acetic acid based, environmentally friendly in normal use**
- **In 10 litre cans**

Cat. No.	Description	€
331701	Liquid neutralizer NLL 10 litres	NC -

Broxal® classic regenerating salt



- **Compacted salt for water softening in dishwashers**
- **Suitable for all types of dishwashers**

Cat. No.	Description	€
155006	Broxal regenerating salt for dishes 2 kg	NC -

Lancer Regenerating Salt



Very high purity white salt, without insolubles that could block the softener. The packaging per kilo corresponds exactly to the quantity needed to fill the salt jar (thus avoiding random dosing). The packaging per unit allows the granulometry to be preserved in order to facilitate its dissolution during a cycle (without risk of agglomeration).

- **Allows easier handling (built-in handle, easier to store) and reduction of packaging (ISO 14001 standard)**

Cat. No.	Description	€
331705B	Regenerating salt - 2 kg	NC -

Regenerating salt for automatic washing machines



Softens the washing water, for a residue-less cleaning.

- **Sodium chloride**
- **Odourless white granule**

Cat. No.	Description	€
711667	Regenerating salt, 5 l bucket	NC -

NEOSEL Regenerating Salt



Recrystallising sodium chloride obtained by evaporation under vacuum and mechanically compressed. Specially developed to obtain a very high purity brine.

Cat. No.	Description	€
148608	NEOSEL Regenerating Salt - 5 kg bucket	NC -