Whatman™-Extraktionspatrone Einfache Dicke



- Für Soxhlet-Extraktionsautomaten
- Hergestellt aus hochwertiger Alphazellulose
- Temperaturbeständig bis 120 °C
- Ausgezeichnete mechanische Festigkeit und Retention

Art. Nr.	Innen-ø x H (mm)	Dicke (mm)	VE	CHF/ VE
036993	10 x 50	1	25	AA -
2800185	18 x 55	1	25	AA -
037228	19 x 90	1	25	AA -
037229	22 x 65	1	25	AA -
036922	22 x 80	1	25	AA -
037058	25 x 100	1	25	AA -
036579	25 x 80	1	25	AA -
036904	26 x 60	1	25	AA -
2800280	28 x 100	1	25	AA -
036948	28 x 120	1	25	AA -
037233	28 x 80	1	25	AA -
037083	30 x 100	1	25	AA -
037167	30 x 80	1	25	AA -
036949	33 x 100	1	25	AA -
28-4189-78	33 x 118	1	25	AA -
036578	33 x 80	1	25	AA -
036974	33 x 94	1	25	AA -
036907	37 x 130	1	25	AA -
036887	41 x 123	1	25	AA -
037191	43 x 123	1	25	AA -
037112	60 x 180	1	25	AA -

Whatman™-Extraktionspatrone Doppelte Dicke Grad 603



- Für Soxhlet-Extraktionsautomaten
- Analyse von organischen Verbindungen
- Hergestellt aus hochwertiger Zellulose
- Ausgezeichnete mechanische Festigkeit und Retention

Art. Nr.	Innen-ø x H (mm)	Dicke (mm)	VE	CHF/ VE
037231	25 x 60	1,5	25	AA -
036909	25 x 80	1,5	25	AA -
036923	25 x 100	1,5	25	AA -
037161	26 x 60	1,5	25	AA -
037143	27 x 80	1,5	25	AA -
037232	28 x 60	1,5	25	AA -
037234	28 x 80	1,5	25	AA -
037038	28 x 100	1,5	25	AA -
036921	30 x 80	1,5	25	AA -
037235	30 x 100	1,5	25	AA -
037237	33 x 60	1,5	25	AA -
036592	33 x 80	1,5	25	AA -
037238	33 x 90	1,5	25	AA -
037057	33 x 94	1,5	25	AA -
037174	33 x 100	1,5	25	AA -
036970	33 x 118	1,5	25	AA -
036842	33 x 130	1,5	25	AA -
036940	33 x 205	1,5	25	AA -
036901	34 x 130	1,5	25	AA -
036938	35 x 150	1,5	25	AA -
037152	41 x 123	2	25	AA -
036939	41 x 145	2	25	AA -
800314	22 x 60	2	25	AA -

Weitere Patronen sind auf Anfrage erhältlich

FILTRATION

Extraktionspatronen

Zur Extraktion von Fett.

- Aus Zellulose
- Dicke 1 mm

Art. Nr.	Außen-Ø x Außen-L. (mm)	VE	CHF/ VE
074185	16 x 100	25	AA -
074186	19 x 90	25	AA -
074187	22 x 60	25	AA -
074168	22 x 80	25	AA -
074188	25 x 80	25	AA -
074180	26 x 60	25	AA -
074189	27 x 80	25	AA -
175028	28 x 80	25	AA -
074190	28 x 100	25	AA -
-	-	-	-
074191	30 x 80	25	AA -
074192	30 x 118	25	AA -
074193	34 x 130	25	AA -
074194	37 x 130	25	AA -
074195	41 x 123	25	AA -
-	-	-	-

Art. Nr.	Außen-Ø x Außen-L. (mm)	VE	CHF VE
074197	46 x 130	25	AA -
074198	48 x 145	25	AA -
074199	50 x 160	25	AA -
074169	33 x 80	25	AA -
175000	33 x 100	25	AA -
175001	41 x 150	25	AA -
-	-	-	-
175003	20 x 80	25	AA -
-	-	-	-
-	-	-	-
175005	33 x 94	25	AA -
-	-	-	-
-	-	-	-
-	-	-	-
-	-	-	-
175010	30 x 100	12	AA -

Art. Nr.	Außen-Ø x Außen-L. (mm)	VE	CHF/ VE
175011	30 x 150	25	AA -
175012	25 x 100	25	AA -
175013	19 x 40	25	AA -
175014	33 x 118	25	AA -
175015	43 x 123	25	AA -
175016	35 x 150	25	AA -
175017	30 x 130	25	AA -
175018	28 x 22	25	AA -
175019	28 x 120	25	AA -
175020	33 x 60	25	AA -
175021	58 x 170	25	AA -
175022	70 x 240	25	AA -
175023	80 x 200	25	AA -
175024	80 x 250	25	AA -
175025	25 x 60	25	AA -
175026	48 x 180	25	AA -
175027	33 x 205	25	AA -

* Innen-Ø

Preisänderung vorbehalten Besuchen Sie www.milian.com