
Magnetische Rührstäbe Spinbar

➤ Aus reinem PTFE, FDA-zugelassen ➤ Autoklavierbar

Achteckige Stäbe mit Drehring






➤ Schwenkring begrenzt die Reibung und verbessert die Effizienz
➤ Achteckige Form für starke Turbulenzen

L x ø (mm)	Farben: Weiß	Rot	Gelb	Blau	CHF
13 x 3	045902	-	-	-	AA -
13 x 8	045907	045903	045904	045906	AA -
13 x 10	046088	-	-	-	AA -
16 x 8	045911	045912	045913	045914	AA -
22 x 8	045915	045916	045917	045918	AA -
25 x 8	045919	045920	045921	045922	AA -
25 x 10	045908	-	-	-	AA -
28 x 8	046091	-	-	-	AA -
38 x 8	045923	045924	045925	045926	AA -
38 x 10	045909	-	-	-	AA -
41 x 8	045931	-	-	-	AA -
51 x 8	045927	045928	045929	045930	AA -
51 x 10	045910	-	-	-	AA -
64 x 8	045932	-	-	-	AA -
64 x 10	045905	-	-	-	AA -
75 x 13	045933	045934	045935	045936	AA -

Mini-Rührstäbchen






➤ Für kleine Gefäße
➤ Seitenstreifen zur Erhöhung der Turbulenzen

L x ø (mm)	Farben: Weiß	Rot	Gelb	Blau	CHF
2 x 2	045945	045947	045946	045948	AA -
3 x 3	045957	045958	045959	045960	AA -
5 x 2	045949	045950	045951	045952	AA -
6 x 3	045961	045962	045963	045964	AA -
7 x 2	045953	045954	045955	045956	AA -
8 x 1,5	045937	045938	045939	045940	AA -
8 x 3	095893	-	-	-	AA -
10 x 3	045965	045966	045967	045968	AA -
12,7 x 3	045969	045970	045971	045972	AA -
15 x 1,5	045941	045943	045942	045944	AA -
20 x 3	095894	-	-	-	AA -
30 x 3	095895	-	-	-	AA -

Kostengünstige Seitenstreifenstäbe







L x ø (mm)	Art. Nr.	CHF
12 x 3	095958	AA -
25 x 8	095959	AA -
40 x 8	095960	AA -
50 x 8	095961	AA -

➤ Seitenstreifen zur Erhöhung der Turbulenzen

Zylindrische Stäbe mit Schwenkring






➤ Schwenkring begrenzt die Reibung und verbessert die Effizienz
➤ Seitenstreifen zur Erhöhung der Turbulenzen

L x ø (mm)	Art. Nr.	CHF	L x ø (mm)	Art. Nr.	CHF
8 x 3	046010	AA -	30 x 6	046073	AA -
12 x 4,5	046065	AA -	35 x 6	046074	AA -
12 x 6	046066	AA -	40 x 8	046076	AA -
15 x 4,5	046068	AA -	45 x 8	046077	AA -
15 x 8	046069	AA -	50 x 8	046078	AA -
20 x 6	046070	AA -	60 x 10	046080	AA -
25 x 6	046071	AA -	70 x 10	046082	AA -
25 x 10	046072	AA -	75 x 13	046083	AA -

Zylindrische Stäbe







➤ Perfekt zylindrisch und glatt
➤ Begrenzte Turbulenzen

L x ø (mm)	Art. Nr.	CHF	L x ø (mm)	Art. Nr.	CHF
6 x 3	046030	AA -	30 x 6	046014	AA -
8 x 3	046031	AA -	30 x 7	046043	AA -
10 x 6	046032	AA -	30 x 8	045994	AA -
12 x 4,5	046033	AA -	35 x 6	046044	AA -
13 x 6	046034	AA -	40 x 7	046045	AA -
13 x 8	045992	AA -	40 x 8	045995	AA -
15 x 4,5	046037	AA -	45 x 8	045996	AA -
15 x 6	046038	AA -	50 x 8	045997	AA -
20 x 6	046039	AA -	57 x 27	045847	AA -
20 x 7	046040	AA -	60 x 10	045845	AA -
20 x 8	046041	AA -	70 x 10	045846	AA -
25 x 5	046042	AA -	80 x 10	045844	AA -
25 x 6	046013	AA -	108 x 27	045848	AA -
25 x 8	045993	AA -	159 x 27	045849	AA -

Stäbe mit Seitenstreifen





➤ Seitenstreifen zur Erhöhung der Turbulenzen

L x ø (mm)	Art. Nr.	CHF	L x ø (mm)	Art. Nr.	CHF
10 x 6	045973	AA -	35 x 7	046055	AA -
12 x 4,5	046046	AA -	35 x 10	045839	AA -
13 x 8	046047	AA -	40 x 7	046056	AA -
15 x 4,5	046048	AA -	40 x 8	045978	AA -
15 x 6	045974	AA -	45 x 8	-	-
15 x 8	046049	AA -	50 x 7	046058	AA -
20 x 6	046012	AA -	50 x 8	045979	AA -
20 x 7	045975	AA -	50 x 10	046059	AA -
25 x 6	046050	AA -	60 x 7	046060	AA -
25 x 8	045976	AA -	60 x 10	045980	AA -
25 x 10	046052	AA -	65 x 13	046062	AA -
30 x 6	046053	AA -	70 x 10	046063	AA -
30 x 7	045977	AA -	75 x 13	046064	AA -
35 x 6	-	-	80 x 7	045981	AA -