

Barreaux d'agitation magnétique Spinbar (suite)

Barreaux forme ovoïde






➤ **Adapté aux ballons et autres récipients à fond rond**

l x Ø (mm)	Pour ballons	Réf.	CHF
10 x 5	-	045982	NC -
15 x 6	-	045983	NC -
19 x 9,5	10 ml à 100 ml	045984	NC -
25 x 12	100 ml à 200 ml	045985	NC -
30 x 16	300 ml à 500 ml	045986	NC -
35 x 16	300 ml à 500 ml	045987	NC -
40 x 13	500 ml à 1000 ml	046097	NC -
40 x 20	500 ml à 1000 ml	045988	NC -
50 x 20	2000 ml à 3000 ml	045989	NC -
64 x 20	3000 ml à 5000 ml	045990	NC -
70 x 20	12 l à 22 l	045991	NC -

Barreaux forme carrée

➤ **Fortes turbulences, particulièrement à basse vitesse**






l x Ø (mm)	Réf.	CHF
12 x 4	095878	NC -
25 x 5,5	095879	NC -
35 x 6	095880	NC -
50 x 7,5	095881	NC -

Barreaux forme profilée






➤ **Action similaire aux formes ovales ou à pivot**

l x Ø (mm)	Réf.	CHF	l x Ø (mm)	Réf.	CHF
10 x 4	095931	NC -	40 x 8	095937	NC -
15 x 5	095932	NC -	45 x 8	095938	NC -
20 x 7	095933	NC -	50 x 8	095939	NC -
25 x 8	095934	NC -	60 x 8	095941	NC -
30 x 8	095935	NC -	70 x 8	095942	NC -
35 x 8	095936	NC -	80 x 10	095943	NC -

Barreaux cylindriques en verre

➤ **Permet d'éviter le contact avec le PTFE**






l x Ø (mm)	Réf.	CHF
6 x 5	095910	NC -
12 x 5	095911	NC -
25 x 6	095912	NC -
45 x 8	095913	NC -
60 x 8	095914	NC -

Barreaux forme triangulaire





➤ **Particulièrement adapté à la dissolution de solides et aux sédiments**


➤ **Le frottement de la base provoque une action de raclement sur le fond**

l x Ø (mm)	Réf.	CHF	l x Ø (mm)	Réf.	CHF
12 x 6	095882	NC -	50 x 12	085887	NC -
20 x 8	095883	NC -	55 x 14	085888	NC -
25 x 8	095884	NC -	80 x 18	085889	NC -
25 x 14	046087	NC -	110 x 36	085890	NC -
35 x 9	095885	NC -	136 x 36	085891	NC -
40 x 14	085886	NC -	-	-	-

Barreaux forme disque

➤ **Adapté aux tubes à essais**






l x Ø (mm)	Réf.	CHF
6 x 9	098889	NC -
6 x 10	095890	NC -
10 x 20	098891	NC -
12 x 30	095892	NC -

Barreaux à pivot amovible





➤ **Utilisable avec ou sans pivot**


➤ **Parfaitement cylindriques et lisses**

l x Ø (mm)	Réf.	CHF	l x Ø (mm)	Réf.	CHF
12 x 8	095919	NC -	51 x 8	095925	NC -
25 x 8	095920	NC -	51 x 10	095926	NC -
32 x 8	095921	NC -	75 x 12	095927	NC -
38 x 8	095922	NC -	102 x 16	095928	NC -
38 x 10	095923	NC -	127 x 16	095929	NC -
42 x 10	095924	NC -	150 x 19	095930	NC -
45 x 8	095975	NC -	-	-	-

Barreaux forme croix

➤ **Très stable, provoque de fortes turbulences**





l x Ø (mm)	Réf.	CHF
10 x 5	095951	NC -
20 x 8	095952	NC -
25 x 9	095953	NC -
30 x 10	095954	NC -
38 x 11	095955	NC -
50 x 15	095956	NC -
60 x 20	095957	NC -

Prix indicatifs. Consultez www.milian.com