Corning® Bar Coded Storage Tubes

A system solution for Automated Sample Management

CORNING

Corning 1D/2D storage line is a superior storage solution designed to provide maximum identification.





Features

- **Maximum Identification** – Synchronized 2D and linear bar code – along with marking spot
- **Maximum Information** – 14 x 14 dot 2D bar code
- **Bar Code Stability** – Laser-etched, not an attached label
- **Compound Compatibility** – Polypropylene construction throughout,making it very inert

	96 1.3 mL 1D/2D Tube	384 Round Tube	384 Square Tube
Tube material	Polypropylene	Polypropylene	Polypropylene
Temperature range	-196°C to +121°C (screw cap) -80°C to +110°C (TPE cap)	-100°C to 110°C (no cap)	-80°C to +110°C (cap)
Tube length	43.9 mm	14.0 mm	14.0 mm
Tube height	52.7 mm with screw cap	16.9 mm with cap	–
	47.5 mm with TPE cap	17.2 mm with tubes, caps, no lid	17.5 mm with tubes, heat seal, no lid
	53.7 mm with tubes, caps, no lid	23.0 mm with lid (no cap)	23.1 mm with tubes, heat seal, and lid
	58.3 mm with tubes, caps, and lid	–	–
Volume	1.3 mL	110 µL	120 µL
Working volume	1.0 mL	70 µL	100 µL
1D codes	Linear, human readable	Code 128 format	–
2D code	14 x 14 dots, 14 digits	14 x 14 dots, 14 digits	14 x 14 dots, 14 digits

Ordering Information

Cat. No.	Description	Qty/Pk	Pk/Cs	Qty/Cs
8500	96 1D/2D bar coded storage tubes, 1.3 mL, without screw caps, bulk pack	96	10	960
8501	96 1D/2D bar coded storage tubes, 1.3 mL, without screw caps, racked with cover	96	10	960
8502	96 1D/2D bar coded storage tubes, 1.3 mL, with screw caps, bulk pack	96	10	960
8503	96 1D/2D bar coded storage tubes, 1.3 mL, with screw caps, racked with cover	96	10	960
8504	96 storage tube screw cap, PP with O-ring, bulk pack	96	50	4,800
8505	96 storage tube screw cap, PP with O-ring, on mat, 96 caps/mat, clear with cover	480	10	4,800
8506	TPE (thermoplastic elastomer) 96 storage tube septum cap, bulk pack	96	100	9,600
8507	TPE 96 storage tube septum cap on mat, 96 caps/mat	960	5	4,800
8508	Empty racks with lids for 96 screw or TPE capped storage tubes	5	2	10
8509	384 2D bar coded storage tubes, round, without plug caps, racked	384	20	7,680
8510	384 2D bar coded storage tubes, square, racked	384	20	7,680
8511	384 storage tube plug cap for round tube, bulk pack	384	100	38,400
8512	384 storage tube plug cap on mats for round tube, 384 caps/mat	3,840	10	38,400
8513	Empty racks with lids for 384 round or square storage tubes	1	20	20

DOMINIQUE DUTSCHER SAS

DOMINIQUE DUTSCHER SAS

How to Purchase: For specific availability in your region and purchasing options, terms and conditions of sale, customer/product support, and certificates, visit **www.corning.com/how-to-buy**.

Warranty/Disclaimer: Unless otherwise specified, all products are for research use or general laboratory use only.* Not intended for use in diagnostic or therapeutic procedures. Not for use in humans. These products are not intended to mitigate the presence of microorganisms on surfaces or in the environment, where such organisms can be deleterious to humans or the environment. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications. *For a listing of US medical devices, regulatory classifications or specific information on claims, visit www.corning.com/resources.

Corning's products are not specifically designed and tested for diagnostic testing. Many Corning products, though not specific for diagnostic testing, can be used in the workflow and preparation of the test at the customers discretion. Customers may use these products to support their claims. We cannot make any claims or statements that our products are approved for diagnostic testing either directly or indirectly. The customer is responsible for any testing, validation, and/or regulatory submissions that may be required to support the safety and efficacy of their intended application.

CORNING

Corning Incorporated
Life Sciences

www.corning.com/lifesciences

NORTH AMERICA
t 800.492.1110
t 978.442.2200

ASIA/PACIFIC
Australia/New Zealand
t 61 427286832
Chinese Mainland
t 86 21 3338 4338
India
t 91 124 4604000
Japan
t 81 3-3586 1996
Korea
t 82 2-796-9500
Singapore
t 65 6572-9740
Taiwan
t 886 2-2716-0338

EUROPE
CSEurope@corning.com
France
t 0800 916 882
Germany
t 0800 101 1153
The Netherlands
t 020 655 79 28
United Kingdom
t 0800 376 8660

All Other European Countries
t +31 (0) 206 59 60 51

LATIN AMERICA
grupoLA@corning.com
Brazil
t 55 (11) 3089-7400
Mexico
t (52-81) 8158-8400

The information contained within is accurate as of the date of publication and subject to change without notice.
For a listing of trademarks, visit www.corning.com/trademarks. All other trademarks are the property of their respective owners.
© 2010, 2025 Corning Incorporated. All rights reserved. 9/25 CLS-MP-054 REV1