


A1


This drawing and/or any technical information relating to it issued in strict confidence and may not be reproduced or used for manufacture without written permission. Upon receipt of which it is the users responsibility to seek confirmation of any dimensions and production tolerances









ALTERNATIVE POSITIONS FOR CUT CORNER





SECTION B-B



VIAL FORM 96pos

NOTE

ALL FACES TO BE HAVE A SMOOTH GLOSS FINISH (HIGHLY POLISHED ON TOOL) UNLESS OTHERWISE STATED

Weight 21g

POLYMIQUE DUTSCHER SAS

AZENTA LIFE SCIENCES www.azenta.com	CUSTOMER	DRAWN ZOLTAN E
DO NOT SCALE	CUSTOMER DRG No / ISSUE	DATE 23 01 2011
IF IN DOUBT ASK	MATERIAL POLYPROPYLENE Purell HP371P	CHECKED PAUL D
	TOLERANCES GENERAL TOLERANCES (UNLESS OTHERWISE STATED) 0 ±0 10; 0 0 ±0 05; 0 0 0 ±0 002; ANGULAR ±0°30' (ALL DIMENSIONS IN MILLIMETRES)	DATE
		SCALE 2 1
		TITLE One piece semi-skirted plate
		DRAWING No Marketing drawing
		SHEET No 1 OF 1 / ISSUE D