

REVISIONS			
REV	ECN	DESCRIPTION	DATE
AA		V-shaped web to be 0.008. (was 0.004) Modify "+" detail to "-" detail	5/5/00





NOTES: UNLESS OTHERWISE SPECIFIED
1. MATERIAL: SILICON RUBBER

Tolerances

0–0.4 in	± 0.004 in
0.4–0.8 in	± 0.005 in
0.8–1.2 in	± 0.008 in
1.2–1.6 in	± 0.01 in
1.6–2.0 in	± 0.012 in
2.0–4.0 in	± 0.024 in
>4.0 in	± 0.04 in

PROPRIETARY AND CONFIDENTIAL

THIS DRAWING IS THE PROPERTY OF AXYGEN SCIENTIFIC INCORPORATED, UNION CITY, CALIFORNIA, AND MAY NOT BE REPRODUCED, PUBLISHED, OR DISCLOSED TO OTHERS WITHOUT WRITTEN CONSENT OF THE COMPANY. IF LOANED IT IS SUBJECT TO RETURN ON DEMAND AND MAY NOT BE USED IN ANY WAY DIRECTLY OR INDIRECTLY DETRIMENTAL TO THE COMPANY.

UNLESS OTHERWISE SPECIFIED
DIMENSIONS ARE IN INCHES
TOLERANCES ARE:
ANGLES: +/-
DECIMAL: +/-
FRACTIONS: +/-

APPROVALS	DATE
DRAWN CRN	08/16/00
CHECKED	
ENGR	

e:\drg\app\sdrg\sl0062aa

AXYGEN SCIENTIFIC
33170 CENTRAL AVENUE UNION CITY, CA 94587 USA

DWG TITLE: AM-96-PCR-RD

SIZE	DWG NO.	REV
B	SL0062	AA

SCALE 1 : 1, 4 : 1 SHEET 1 OF 1

DOMINIQUE DUTSCHER SAS