Nalge Nunc International Corporation,
a Part of **Thermo Fisher Scientific Inc.**
75 Panorama Creek Drive
Rochester, NY 14625
1 (585) 586-8800
www.thermofisher.com

March 6, 2025

Subject: Statement on Animal Derived Components for Thermo Scientific™ Nalgene™ Catalog Number Series 4103-XXXX and 4110-XXXX

Dear Customer:

Per your request, below is the Thermo Fisher Scientific statement on Animal Derived Components for Nalgene catalog number series 4103-XXXX and 4110-XXXX. Thermo Fisher Scientific manufacturing uses one raw material in the production of the flasks.

The flasks are manufactured using polycarbonate (PC) resin, part number 8-0056-53. The PC resin supplier has indicated that one of the additives used in the formulation of the resin is a derivative of fatty acid, made from beef tallow. However, the tallow is rendered inactive in accordance with the European Union Directive 76/768/EEC, under 200oC/40 bar/20 min condition.

Thermo Fisher Scientific does not reserve production equipment for raw materials that are free of animal origin ingredients. All resins and colorants used at the manufacturing site comply with the latest revision of EMA/410/01 section 6.4.

Sincerely,



Jessica Garigen
Regulatory Specialist
Thermo Fisher Scientific
Labware and Specialty Plastics