
PRODUCT DRAWING



ALUMINUM FOIL WITH ACRYLIC ADHESIVE ON A POLY-COATED PAPER LINER

BACKSCORES THROUGH LINER ONLY

PERFORATIONS THROUGH FILM ONLY PERF: START CUT .060"(1.52) / TIE .032"(0.81) / CUT .050"(1.27)

DOMINIQUE DUTSCHER SAS

UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES. "( )" DIMENSIONS ARE IN MILLIMETERS. UNSPECIFIED TOLERANCES ARE ± 0.015.						
	T.CARROLL	DATE: 5/19/14	DESC: AlumaSeal II™, AF-NL-100			
CHECKED	B.CAHILLY	DATE: 5/21/14				
APPROVED	B.CAHILLY	DATE: 5/21/14	SCALE: 1 : 1	DWG #: MFG-DWG-405	REV: A	SHEET: 1 OF 1