


A1

This drawing and/or any technical information relating to it issued in strict confidence and may not be reproduced or used for manufacture without written permission. Upon receipt of which it is the users responsibility to seek confirmation of any dimensions and production tolerances.




CAP DETAIL 10:1





DOMINIQUE DUTSCHER SAS

NOTE:
ALL FACES TO HAVE A SMOOTH GLOSS FINISH (HIGHLY POLISHED ON TOOL)



ACTUAL SIZE

ISSUE	DATE	AMENDMENT

AZENTA LIFE SCIENCES
www.azenta.com

DO NOT SCALE
IF IN DOUBT ASK

CUSTOMER:		DRAWN:	ZE.
CUSTOMER DRG. No.:	ISSUE:	DATE:	08.01.2020
MATERIAL:	POLYPROPYLENE	CHECKED & SIGNATURE:	PP.
TOLERANCES: GENERAL TOLERANCES (UNLESS OTHERWISE STATED)	0 ±0.15; 0.0 ±0.05; 0.00 ±0.02; ANGULAR ±0°30' (ALL DIMENSIONS IN MILLIMETRES)	DATE:	
		SCALE:	5:1

TITLE: Strip of 8 domed optical cap Marketing drawing
DRAWING No.: 4ti-0752
SHEET No.: 1 OF 1
ISSUE: D