


A0

This drawing and/or any technical information relating to it issued in strict confidence and may not be reproduced or used for manufacture without written permission. Upon receipt of which it is the users responsibility to seek confirmation of any dimensions and production tolerances.







VIAL FORM 96pos
SCALE 10 1



SECTION A-A





4Ti-0900/C clear frame clear well
4Ti-0901 clear frame white well

DOMINIQUE DUTSCHER SAS

ISSUE	DATE	AMENDMENT

AZENTA LIFE SCIENCES
www.azenta.com

DO NOT SCALE
IF IN DOUBT ASK

CUSTOMER		DRAWN	ZOLTAN E
CUSTOMER DRG No	ISSUE	DATE	08 04 2009
MATERIAL		CHECKED	PAUL D
TOLERANCES	0 ±0 25 0 0 ±0 20 0 00 ±0 10 ANGULAR ±0°30'	DATE	
GENERAL TOLERANCES (UNLESS OTHERWISE STATED) (ALL DIMENSIONS IN MILLIMETRES)		SCALE	

TITLE: FrameStar Semi Skirted CBL
DRAWING No: 4Ti-0900
SHEET No: 1 OF 1
ISSUE: A