Date: 09/19/2023

From: Whirl-Pak®

Subject: Endotoxins in Whirl-Pak® Bags

Because endotoxins can come from anywhere at any time, Whirl-Pak® bags are not tested for their presence. The bags are sterilized after manufacturing with ethylene oxide gas which will eliminate contamination from endotoxins, and protect the sample once it is placed inside the bag. Parametric Release Process is used to verify sterility and documentation is available anytime from the whirl-pak website at www.whirl-pak.com

If you have any questions, or need assistance, please email our Customer Support team at contact@whirl-pak.com. For Frequently asked questions please visit our website www.whirl-pak.com/frequently-asked-questions