• Diversified Biotech •

Product Specification Sheet

Product Name: Tough-Spots® on Rolls – Label part numbers beginning with T-SPOTS and T-SPOTS-50 and/or SPRL, BHTS and , Cryo-Tags® and Cryo-Babies® (on rolls)

Face stock : 3.5 mils (100.0 microns) Flexible White Vinyl Film

Adhesive : Acrylic – 1.0 mils (25.0 microns) #300 "Hi-Strength" Clear Acrylic Adhesive

Liner : 3.2 mils (77.5 microns) 50# White Semi-Bleached Densified Kraft Paper Liner*

Temperature Resistance: -196°C to 80°C

Shelf Life: One year from date of manufacture. Store at 72°F and 50% relative humidity.
HARM Code - 3919.10.1000
Country of Origin- USA
USMCA level B
EAR99

SDS (US OSHA)
This product is an article and therefore is not subject to the requirements of the US Occupational Safety and Health
Administration's (OSHA) Hazardous Communications Standard 29 CFR 1910.1200(b)(6)(v) to provide a Safety Data Sheet
(SDS).

China RoHS
This product or part does not exceed the maximum concentration values (MCVs) in China RoHS GB/T 26572-2011
(Requirements of Concentration limits for certain restricted substances in electrical and electronic products), which means
that each of the homogeneous materials within this product is below the following MCVs: (a) 0.1% (by weight) for lead and
its compounds, mercury and its compounds, hexavalent chromium and its compounds, polybrominated biphenyls or
polybrominated diphenyl ethers; and (b) 0.01% (by weight) for cadmium and its compounds.

CONEG
This product complies with the US CONEG (Coalition of Northeastern Governors) Model Toxics in Packaging Legislation as
there is no intentionally added lead, cadmium, hexavalent chromium or mercury, and the sum total concentration of these
substances in the product does not exceed 100 ppm by weight.

Conflict Minerals
Conflict Minerals, which the U.S. Securities and Exchange Commission ("SEC") has defined as gold, columbite-tantalite

DOMINIQUE DUTSCHER SAS

(coltan), cassiterite, wolframite, or their derivatives (tin, tantalum, or tungsten), are not contained in or are not “necessary to
the functionality or necessary to the production” of the above-listed product, as the term “necessary to the functionality or the
production” is defined under the SEC’s Conflict Minerals Rule. 77 Fed. Reg. 56274 (Sept. 12, 2012).

EU Azo Colourant Directive

This product does not contain azo dyes which, by reductive cleavage of one or more azo groups, may release one or more of
the aromatic amines listed in the Appendix I, in detectable concentrations above 30 ppm in the finished article or in the dyes.
parts thereof, according to the testing method established in accordance with Article 2a of this European Commission
Directive 2002/61/EC (2004/21/EC).

EU CMRs

This product does not contain an ingredient at >= 0.1% that is classified as a Category 1A or 1B carcinogen, mutagen, or
reproductive toxicant according to Annex VI of European Regulation 1272/2008 on Classification, Labelling, and Packaging
of Dangerous Substances and Mixtures or that is classified by 3M or its vendors as a Category 1A or 1B carcinogen, mutagen,
or reproductive toxicant according to the criteria laid out in European Regulation 1272/2008 on Classification, Labelling, and
Packaging of Dangerous Substances and Mixtures.

EU REACH

This product is an article, without intended release of a chemical substance, under the Regulation No 1907/2006 of the
European Parliament and the Council concerning the Registration, Evaluation, Authorisation and Restriction of Chemicals
(REACH) (refer to REACH, Article 3(3)). It is not a chemical preparation. Therefore, it is not subject to the (pre)-registration.
or the registration process. It does not require a safety data sheet.

EU REACH

This product, including any article that the product is composed of, does not contain at greater than 0.1% by weight a
Substance of Very High Concern (SVHC) substance identified according to Article 59 of REACH. This declaration reflects.
the substances on the candidate SVHC list, effective June 2023.

EU RoHS

This product does not exceed the maximum concentration values (MCVs) set under EU Directive 2011/65/EU (RoHS
recast/RoHS 2), as stated in Annex II to that directive. This means that each of the homogenous materials within this product.
does not exceed the following MCVs: (a) 0.1% (by weight) for lead, mercury, hexavalent chromium, polybrominated
biphenyls or polybrominated diphenyl ethers; and (b) 0.01% (by weight) for cadmium.

EU RoHS Phthalates

This product does not exceed the maximum concentration values (MCVs) for phthalates set under EU Directive 2011/65/EU
(RoHS recast/RoHS 2), as amended by EU 2015/863, which applies to finished EEE after July 22, 2019 for Category 1-7, 10-

DOMINIQUE DUTSCHER SAS

11 products and after July 22, 2021 for Category 8 and 9 products. This means that each of the homogeneous materials within
this product does not exceed the MCV of 0.1% (by weight) for each of the following phthalates: DEHP, BBP, DBP, and
DIBP.

Fluorinated Greenhouse Gas
A Fluorinated Greenhouse Gas (or F-Gas) listed under Annex I and II of Regulation (EU) No 517/2014 of the European
Parliament and of the Council of 16 April 2014 is not intentionally added to the product.

GADSL
This product does not intentionally contain a Declarable or Prohibited substance(s) at or above the reportable threshold levels
in the Global Automotive Declarable Substance List (GADSL). This declaration reflects the substances on the list as of the
date of this certificate. *7604FP: 3M(tm) White Vinyl Label Material is available in the public domain of IMDS as a Semi-
Component Module entry with ID 424255226.*

Halal
This product has not been certified Halal.

Kosher
This product has not been certified Kosher.

Stockholm Agreement
This product does not contain Persistent Organic Pollutants (POPs) at or above applicable thresholds per Annexes A, B, and
C of the Stockholm Convention, May 2004 and subsequent amendments.

Sustainability Advantage: Recycled content
This product does not contain recycled content.

TSCA Section 6
This product is not known to contain 2,4,6-Tri-tert-butylphenol (CAS 732-26-3).

TSCA Section 6
This product is not known to contain Decabromodiphenyl Ether (Deca-BDE) (CAS 1163-19-5).

TSCA Section 6
This product is not known to contain Hexachlorobutadiene (HCBD) (CAS 87-68-3).

TSCA Section 6
This product is not known to contain Pentachlorothiophenol (PCTP) (CAS 133-49-3).

TSCA Section 6
This product is not known to contain Phenol, isopropylated phosphate (3:1) (PIP (3:1)) (CAS 68937-41-7).

Chemicals and/or Compounds of Interest

Asbestos : Not known to contain.
Halogenated Compounds : Contains.
Materials of Human or Animal Origin : Not known to contain.

Per- and Polyfluoroalkyl Substances (PFAS)

PFAS refers to a broad range of thousands of materials with varying properties and characteristics, which may include high
resistance to oil, water, temperature, chemicals, and fire, which makes certain PFAS critical to the functioning of many
modern products.

Per- and polyfluoroalkyl substances (PFAS) : Not intentionally added.
Perfluorohexane sulfonate (PFHxS) : Not intentionally added.

DOMINIQUE DUTSCHER SAS

Perfluorohexanoic acid (PFHxA) (CAS 307-24-4) : Not intentionally added.
Perfluorooctanesulfonic Acid (PFOS) (CAS 1763-23-1) : Not intentionally added.
Perfluorooctanoic Acid (PFOA) (CAS 335-67-1) : Not intentionally added.

Diversified Biotech Inc. recommends the end user evaluate material for specific requirements. 08/18//2023

DOMINIQUE DUTSCHER SAS