MultiScreen® Plasmid$_{384}$ Miniprep Clearing Plate

For use in lysate clearing during plasmid miniprep

LSKP 384 50 (50/pk)

Recommended Use

The MultiScreen Plasmid$_{384}$ Miniprep Plate is designed for the clearing of plasmid lysates prior to alcohol precipitation. The filtration step is driven by centrifugation. The product is intended for research use only.

Product Specifications

Basic function	Clearing of bacterial cell lysates	**Filtrate volume**	Up to 100 µL
Mode of operation	Centrifugal	**Solvent compatibility**	Aqueous
Maximum g force	Up to 4000 × g	**Aqueous liquid incubation**	Up to 15 minutes
Well volume	Up to 100 µL	**Collection plate**	SBS standard, 384 well plate

Technical Assistance

For more information, contact the Millipore office nearest you. In the U.S., call **1-800-MILLIPORE** (1-800-645-5476). Outside the U.S., see your Millipore catalogue for the phone number of the office nearest you or go to our web site at www.millipore.com/offices for up-to-date worldwide contact information. You can also visit the tech service page on our web site at www.millipore.com/techservice.

Standard Warranty

Millipore Corporation ("Millipore") warrants its products will meet their applicable published specifications when used in accordance with their applicable instructions for a period of one year from the date of shipment of the products. **MILLIPORE MAKES NO OTHER WARRANTY, EXPRESSED OR IMPLIED. THERE IS NO WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.** The warranty provided herein and the data, specifications and descriptions of Millipore products appearing in Millipore's published catalogues and product literature may not be altered except by express written agreement signed by an officer of Millipore. Representations, oral or written, which are inconsistent with this warranty or such publications are not authorized and if given, should not be relied upon.

In the event of a breach of the foregoing warranty, Millipore's sole obligation shall be to repair or replace, at its option, the applicable product or part thereof, provided the customer notifies Millipore promptly of any such breach. If after exercising reasonable efforts, Millipore is unable to repair or replace the product or part, then Millipore shall refund to the customer all monies paid for such applicable product or part. **Millipore shall not be liable for consequential, incidental, special or any other indirect damages resulting from economic loss or property damages sustained by any customer from the use of its products.**

DOMINIQUE DUTSCHER SAS

©2006. Millipore Corporation. All rights reserved.
Millipore and MultiScreen are registered trademarks of the Millipore Corporation.
PR02214, Rev. A, 09/06

MILLIPORE