MTC Bio Incorporated

Certificate of Quality

Catalog No: **C2000-OB**

Lot No: **24/05/15**

Product Description:	1.5mL Microcentrifuge Tube, sterile, with snap-cap, black	**Material:**	Polypropylene (PP)
Packaging:	500/pk	**Country of Origin:**	Manufactured in China

Quality System Compliance:

This product is manufactured under the ISO 9001:2008 standard. Stringent quality control protocols are followed at every process of production ranging from qualitative raw materials and process of production through final products.

Quality Control Testing:

Products are inspected and controlled through the entire production process in accordance with current applicable product specifications and QC SOP. Inspection records are reviewed and signed by qualified personnel for product release.

Product Specifications:

Measurement:	Length: 39.3mm Diameter: 10.4mm
Performance:	Leak-proof, RCF 28,000xg with appropriate adapter
Sterilization:	This product has been sterilized under ISO 11137:2006 guidelines using Electron beam irradiation. Dose Range: 15.0 - 24.0 kGy
DNase & RNase Free	This product has been tested and found free of any detectable DNase and/or RNase contamination
Non-Pyrogenic:	The acceptance level for this product is 0.25 EU/ml (TAL Gel Clot Method)
Cytotoxicity:	Raw materials have been tested for cytotoxicity and have been shown to be non-toxic. (ISO 10993 testing standard)
BSE/TSE Statement:	This product contains no ingredients of animal origin and no material derived from, or exposed to animals affected by or under quarantine for Transmitting, Animal Spongiform Encephalopathy (TSE)/Bovine Spongiform Encephalopathy (BSE)
Latex Statement	This product does not contain latex
RoHS Compliance:	Raw materials used are in conformity with CONEG standards for heavy metals.
REACH (SVHC) Statement:	Raw materials used are in conformity to REACH standards.

Production Information:

Production date: 15, May 2024
Expiration date: 15, June 2029
Sterility expiration: 15, June 2029

QA Dept. Director: _____

DOMINIQUE DUTSCHER SAS