pipetman®

THE PIPETTING STANDARD

RECOGNIZED RELIABILITY AND PROVEN DURABILITY COMBINED WITH EASY PIPETTING AND HANDLING



Legendary Performance

Accurate and precise, PIPETMAN® air-displacement pipettes are recognized as one of the world standard. They come in a large range of volumes, including single and multichannel models, allowing you to handle various sample types with confidence.

Easy Pipetting

PIPETMAN feels like PIPETMAN® Classic but with lower forces for enhanced pipetting habits with a standard pipette.

Built-to-Last and Convenient

Made of PVDF and stainless steel, PIPETMAN offers premium durability. Built-to-last, it is the economical solution for years of pipetting.

Simple to Maintain

Routine cleaning and maintenance are all that is required to keep PIPETMAN in top condition for years. The PIPETMAN design allows for easy access to user-serviceable parts, along with a fully autoclavable tip holder and ejector.

DOMINIQUE DUTSCHER SAS

GILSON®


Choose Your Pipette Stand

Different pipette stands and holders allow you to store pipettes vertically, protecting them from contamination.

CARROUSEL™



TRIO™



Universal Multichannel Stand




SINGLE™



Pipette Stands	
Description	**Part Number**
CARROUSEL™ Pipette Stand (7 pipettes)	**F161401**
TRIO™ Pipette Stand (3 pipettes)	**F161405**
Universal Multichannel Stand	**F161417**
SINGLE™ Pipette Holder	**F161406**

Choose Your Pipette Kits

PIPETMAN single channel models are available in kits, offering the essential pipetting tools to start an experiment right away.



PIPETMAN® Kits		
Model	**Part Number**	**Models**
PIPETMAN® 4-Pipette Kit	**F167360**	P2 P20 P200 P1000
PIPETMAN® Starter Kit	**F167900**	P20 P200 P1000
PIPETMAN® Micro-volume Kit	**F167800**	P2 P10 P100
PIPETMAN® Large-volume Kit	**F167920**	P5000 P10mL

PIPETMAN® Single Channel				
Model	**PIPETMAN® DIAMOND Tips**		**Part Number**	**Volume**
P2	D10 DL10	DF10ST DFL10ST	**F144054M**	0.2-2 µL
P10	D10 DL10	DF10ST DFL10ST	**F144055M**	1-10 µL
P20	D200	DF30ST	**F144056M**	2-20 µL
P100	D200	DF100ST	**F144057M**	10-100 µL
P200	D200 D300	DF200ST DF300ST	**F144058M**	20-200 µL
P1000	D1000 D1200	DF1000ST DF1200ST	**F144059M**	100-1000 µL
P5000*	D5000		**F144066**	500-5000 µL
P10mL*	D10mL		**F144067**	1-10 mL
PIPETMAN® Multichannel				
P8x10 **P12x10**	D10 DL10	DF10ST DFL10ST	**F144068** **F144069**	1-10 µL
P8x20 **P12x20**	DL10 D200	DFL10ST DF30ST	**F144070** **F144071**	2-20 µL
P8x200 **P12x200**	D200 D300	DF200ST DF300ST	**F144072** **F144073**	20-200 µL
P8x300 **P12x300**	D200 D300	DF200ST DF300ST	**F144074** **F144075**	30-300 µL

Models equipped with a tip ejector are only sold with a stainless steel tip ejector. Plastic tip ejector is available as an accessory. Refer to *PIPETMAN® User's Guide LT801122.*

*P5000 and P10mL do not have tip ejectors.

Gilson maximum permissible errors are guaranteed only when PIPETMAN® pipettes are used with the recommended PIPETMAN DIAMOND Tips.



PIPETMAN DIAMOND Tips are designed to fit your PIPETMAN perfectly and guarantee maximum accuracy and precision.

For ordering information, visit www.gilson.com

DOMINIQUE DUTSCHER SAS

gilson.com/contactus

LT800591/F | ©2023 Gilson S.A.S.
All rights reserved. All trademarks are the property of Gilson, Inc. and its subsidiaries.

GILSON®