WARING COMMERCIAL


50-250 ml Mini Container

# MC3


***Units Marked for Deletion***

**MCC3** *(on 4/2/2002)*
mini cool container


**MMGC3** *(on 4/5/2002)*
medical grade container




**MCC3**
**500139** Blending Assembly




**MC3**
**500139**
Blending Assembly




**MMGC3**
**501380** Blending Assembly
(no longer available)

DOMINIQUE DUTSCHER SAS

| Catalog# | Model |
|---|---|
| MC3 | 5019 |
| *Units Marked for Deletion* | |
| MCC3 | 31CA75 |
| MMGC3 | 31CA85 |

Waring Factory Service Center
314 Ella T. Grasso Ave.
Torrington, CT 06790
Tele. 1-800-269-6640
Fax 860-496-9017
www.waringcommercialproducts.com

| Illustration # | Part # | Description |
|---|---|---|
| 1 | 019448 | Lid |
| 2 | 027862 | Cap Nut *(on units after 7/2/2021)* |
| | 003504 | Cap Nut |
| 3 | 002352 | Blade |
| 4 | 003609 | Stainless Steel Washer |
| 5 | 003508 | Fairprene Washer |
| | 005321 | Small Special Washer (for 501701 Special Blending Assembly) |
| 6 | 010012 | Bearing Holder |
| 7 | 009724 | Nylon Washer |
| 8 | 011196 | Container |
| 9 | 003537 | Stainless Steel Washer |
| 10 | 003509 | Washer /neoprene |
| | 005320 | Large Special Washer (for 501701 Special Blending Assembly) |
| 11 | 015055 | Container Adapter |
| 12 | 009780 | Hex Nut |
| 13 | 003510 | Washer |
| 14 | 002406 | Socket Head Shaft |
| 15 | 500138 | Container and Jacket Assembly (obsolete) |
| 16 | 009762 | Gasket |
| 17 | 005552 | Cap Nut |
| 18 | 016267 | Blade (no longer available) |
| 19 | 015880 | Bearing Holder (no longer available) |
| 20 | 016751 | Hand Nut (obsolete) |
| 21 | 015639 | Socket Head Shaft (obsolete) |

**501701**
MC3 Special Blending Assembly

DOMINIQUE DUTSCHER SAS